PETERS SURGICAL
Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT
Tél. +33 (0)1 48 10 62 62 - Fax +33 (0)1 48 91 22 99
E-mail : peters@peters-surgical.com - www.peters-surgical.com FRANCE

IC5330 SC5330 C5470 F5330 MD5330 H5330
IK5330 SK5330 A5330 G5330 ECS5330 H5330L

Français

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VECTEC®

CE0459

INSTRUMENTS LAPAROSCOPIQUES A USAGE UNIQUE

INCUT / INCUT SMART

CISEAUX MONOPOLAIRES

ENDOGRIP

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STOCKAGE / MANIPULATION / ELIMINATION

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ENDOHOOK
MONOPOLAR HOOK

CAUTION

Please read these Instructions carefully prior to using any single-use laparoscopic instruments. Failure to comply with the following recommendations may entail serious consequences. Laparoscopic surgery should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved in using such surgical technique; notably risks for the patient when manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure. Patients, having a right to information according to Law, should be adequately briefed on the surgical procedure.

These Instructions solely contain information on single-use laparoscopic instruments and their means of use. On no account should they be considered as an instructions' manual for using laparoscopic instruments in general.

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

Rx ONLY (for USA): Caution: US Federal Law restricts this device to sale by or on the order of a physician.

DESCRIPTION

Single-use laparoscopic instruments are fitted with an insulated body of 5 mm in diameter and of varying lengths depending on the model (330 or 470 mm). They are equipped with a monopolar connector of 4 mm in diameter to be connected with a surgical monopolar cable and a high-frequency (HF) generator, expect for forceps without connector ref. ECS5330. They are designed for insertion purposes and to be used with trocars of same diameter or perhaps wider (insofar as equipped with a reducer).

INDICATIONS OF USE

Single-use laparoscopic instruments are solely designed for laparoscopic surgery.

INCUT and INCUT SMART Monopolar scissors are designed for mechanical incision and dissection, monopolar incision and monopolar coagulation via high-frequency current. **ENDOGRIP** Monopolar dissector are designed for mechanical dissection, atraumatic prehension of tissues and monopolar coagulation via high-frequency current. **ENDOGRIP** grasping forceps are designed for mechanical dissection, traumatic or atraumatic prehension of tissues and monopolar coagulation via high-frequency current (expect for ref. ECS5330). **ENDOHOOK** are designed for mechanical dissection, monopolar incision and monopolar coagulation via high-frequency current.

CONTRAINDICATIONS

Use of the single-use laparoscopic instruments is strictly forbidden whenever a contraindication exists for laparoscopic surgery.

Not intended for use on new-born or infant (15 kg min.).

WARNINGS AND PRECAUTIONS

Single-use laparoscopic instruments are designed for single-use only. On no account may they be used a second time over, whatever the reason or the means.

Such instruments should only be handled by qualified practitioners having gained due experience thereof.

All instruments are delivered in sterile packaging. They should never be resterilized, regardless of whether they have been used or not.

Any reuse or reesterilization of this instrument shall affect its original sterility and may lead to infection or to cross-contamination in patients. Reuse or reesterilization may also compromise the integrity and functioning of the instrument, possibly causing injury, sickness, or even death of a patient.

Never use the sterile packaging and its contents in the following cases (whereupon everything should be destroyed):

- Packaging severely crushed or clearly damaged;
- Visible damp marks on the cardboard packaging or on the protective sachet.
- Perforation of the protective sachet.
- Expiry of the use-by date printed on the packaging.
- Damaged insulation sheath.

Damaged instruments should not be used.

WARNING AND PRECAUTIONS FOR ELECTROSURGERY

Electrosurgery is liable to cause direct or indirect scalding. Hence, all necessary precautions should be taken to avoid/reduce any such risks.

Maximum rated voltage: **4,000 at peak voltage** (the generator should be set at a lower voltage).

Ensure compatibility with all other associated instruments; in particular, the connections and voltage supplied by the HF generator should be lower than the rated voltages of the accessories being used in combination (refer to IEC 60601-2-2). Check the instruction for use of the cable and of the HF generator regarding its operating, as well as its power and mode settings.

Only surgical monopolar cables can be connected to PETERS SURGICAL monopolar instruments.

The HF function **should, as far as possible, be used under minimal power and under the constant eye of the user: risk of scalding.**

Electrosurgery may induce electromagnetic interference.

Risk of burning and explosion in the presence of inflammable substances; possible formation of sparks. To avoid the risk of fire, keep electrosurgical instruments away from flammable or oxidizing gases (e.g. O₂, N₂O, endogenous gases, and certain anesthetics), flammable solutions (e.g. certain cleaning and disinfecting agents), and flammable materials (e.g. gauze, cotton and surgical drapes). Use electrosurgical instruments with caution in an oxygen rich environment.

Always make sure that the instruments or cables are not placed directly on the patient, or in a manner to create direct or indirect contact (in particular, when not in use).

Avoid skin-to-skin contact during use of electrosurgical instruments. Current can be transferred and may lead to patient injury or burns.

Electrosurgery may induce muscular contractions.

Immediately after a long use of electrocoagulation, temperature of the jaws may increase enough to damage tissue on contact. Always keep laparoscopic sight of the distal part of the instrument.

Depending on the nature of the tissue and its surface hygrometry, take into consideration the risk of thermal diffusion (risk of tissue necrosis, post-operative bleeding).

Excessive tightening with the jaws can damage the tissues, chose the appropriate type of jaw and use the minimum force needed for manipulation.

Electrical diffusion: in order to prevent damage to tissue not retained in the pincers, do not activate monopolar energy:

- If the tip of the jaws and the surgical surface are not visible,
- If the stem of the jaws comes into contact with the skin tissue,
- In the presence of excessive conductor liquids.

Risks of interaction for patients carrying a pacemaker (fibrillation) or metal implants (scalding).

Monopolar Scissors: do not apply high-frequency (HF) current with the blades open.

INSTRUCTIONS FOR USE

1. Open the multi-pack box and take out the medical device.
2. Make sure that the sterile packaging protecting the medical device is neither open nor damaged. In the slightest doubt, destroy the sterile packaging, together with its contents.
3. Check the use-by date printed on the packaging.
4. Open the sterile packaging according to proper practices.
5. Remove the protective caps and check whether the instrument is in good condition, notably the surface of the isolating sheath. Check the opening/closing of the blades/clamps.
6. Place the instrument on a sterile zone paying attention not to damage the active part through knocking or otherwise.
7. Please take the necessary precautions when inserting or withdrawing the instrument into/from the trocar: risk of damaging the seals of the trocar and of loss of pneumoperitoneum, notably when withdrawing hooks/clamps/dissectors.

Ensure that insulation is intact. Do not use instrument if insulation is damaged or cracked.

Upon completion of surgery, as of when the instrument is no longer inside the trocar, the protective cap should be replaced on the distal part.

STORAGE CONDITIONS / MANIPULATION / DISPOSAL

Storage conditions : Store the unopened sterile package in a cool dry area, with good ventilation away from areas of extreme temperature and humidity.

Transportation conditions: For transportation, the above points should be observed.

Use conditions: Instrument should be used within the following conditions: Temperature between 10°C and 40°C, humidity (RH%) between 30 and 70% and atmospheric pressure conditions between 700 and 1060hPa.

Disposal: After surgery, the instrument should be disposed of in the container reserved for sharp, traumatizing or perforating objects. It is disposed of as per the local or national hospital waste disposal management procedures.

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