INSTRUCTIONS FOR USE: Uterine ElevatOR PRO with OccludOR Balloon™

INDICATIONS FOR USE / CLINICAL BENEFITS
The Uterine ElevatOR PRO with OccludOR Balloon™ is indicated for manipulation of the uterus, and injection of fluids or gasses during laparoscopic procedures including total laparoscopic hysterectomy, laparoscopic assisted vaginal hysterectomy, laparoscopic tubal occlusion, and diagnostic laparoscopy. The Uterine ElevatOR PRO with OccludOR Balloon™ maintains pneumoperitoneum during laparoscopic procedures by sealing the vagina once colpotomy is performed. The device mobilizes the uterus, enhancing visualisation of key structures and anatomical landmarks. In addition, its ability to articulate and rotate the uterus and cervix facilitates the safer and more optimal use of other surgical devices.

CONTRAINdications

Device should not be used in:

- patients with a known or suspected pregnancy
- patients with known or suspected uterine or tubal infection
- patients with a history of pelvic adhesions
- diseases associated with reduced fertility (ART) procedure related to in vitro fertilization
- patients with early stage cervical malignancy undergoing a laparoscopic approach

This device is not intended for use as a Laser Backstop.

CAUTIONS

- The device is intended for single use only. The design of this device may not perform as intended by the manufacturer if it is reused. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device.
- Caution should be exercised in patients with cervical or uterine pathology.
- To prevent air leak, test the intrauterine balloon for proper inflation and deflation of the intrauterine balloon with air.
- Do not insert or remove the Uterine ElevatOR PRO with OccludOR Balloon™ with intrauterine balloon inflated.
- Lubricate the tip and handle of the device before use.
- Check to ensure the correct insertion of the Uterine ElevatOR PRO with OccludOR Balloon™. The cervical os should be at least 5mm (0.2in) in diameter to facilitate safe insertion and to reduce the risk of uterine perforation. Injury to the cervix and/or intrauterine balloon rupture.
- When injecting any liquid media, closely follow the manufacturer's directions for use accompanying the product.
- Use of the Uterine ElevatOR PRO with OccludOR Balloon™ in hysterosalpingography and Rubin's Test should be completed under anesthesia.
- Caution must be used when using energy devices with this device. Any prolonged direct contact may damage devices.
- The device is designed to accommodate no more than 12cc of inflation saline.

WARNINGS - GENERAL

All users must read and be familiar with all instructions, warnings, contraindications and precautions in this package insert before using this device.

- This device should be used only by surgeons trained in proper techniques for intrauterine surgery, laparoscopic surgery, diagnostic pelvic anatomy, and placement of intrauterine retracting instruments.
- As with all intrauterine devices, inspection of the Uterine ElevatOR PRO with OccludOR Balloon™ can result in uterine perforation and subsequent bleeding. The surgeon should examine the patient carefully before and after insertion to ensure proper placement and to prevent uterine perforation or subsequent bleeding and take appropriate clinical steps.

To avoid potential injury to the uterine wall or inadvertent expulsion of the device from the uterus during the procedure, DO NOT undilate the intrauterine balloon. Inflation with 7-10cc of air is recommended to compensate for injected air pressure in the pilot balloon and inflation tube. DO NOT exceed 10cc of air.

DO NOT use any liquids to inflate the intrauterine balloon. Use of fluid could potentially over-stretch and rupture the balloon.

DO NOT inject any dyes / marking fluids directly into the pilot balloon. Dyes / marking fluids should only be injected through the inject port on the rear of the handle. DO NOT attach other devices or lines to the pilot balloon.

DO NOT inject dyes / marking fluids rapidly.

DO NOT attach other devices or lines other than a fuer lock syringe to inject dyes / marking fluid to the injector port.

The Uterine ElevatOR PRO with OccludOR Balloon™ must be inserted along the direction of the uterine cavity, (anteriorly in an antverted uterus, posteriorly in a retroverted uterus) to reduce the risk of uterine perforation. Sound the uterine cavity to assess the uterine depth. Extra caution should be exercised in the case of a very small uterus.

DO NOT use the Uterine ElevatOR PRO with OccludOR Balloon™ as a uterine sound.

Ensure cervical cup is not under tension. Under tensioned cup may cause misplacement of cervical cup around cervical os.

The Uterine ElevatOR PRO with OccludOR Balloon™ is not recommended for use in a large, post-partum uterus as movement of the device may result in uterine wall lacerations and subsequent bleed.

Check the pilot balloon frequently to ensure inflation of the intrauterine balloon. If the intrauterine balloon ruptures, the pilot balloon will not feel firm when squeezed. If the intrauterine balloon has ruptured, stop all manipulation immediately. Remove and replace the Uterine ElevatOR PRO with OccludOR Balloon™ with a new unit.

When in use, the locking assembly must remain locked at all times to ensure the safe positioning of the cervical cup and vaginal cup. If the locking assembly becomes inadvertently open or loose, immediately stop all manipulation and replace the Uterine ElevatOR PRO with OccludOR Balloon™ and dissecting instruments must be avoided during uterine dissection and excision to reduce the risk of patient injury. Damage to the device and/or intrauterine balloon may also occur.

If bending of the Uterine ElevatOR PRO with OccludOR Balloon™ is deemed necessary for use with robotic instruments, all device systems must first be tested prior to use to ensure functionality.

If the surgeon has used the Uterine ElevatOR PRO with OccludOR Balloon™ and the uterus has not been removed, the surgeon should examine the patient carefully for any indications of uterine perforation or uterine wall lacerations and subsequent bleed.

This product contains phthalates that may leach out into the patient during use. Animal studies have shown that high-level exposure to DEHP may affect male fertility and reproductive development.

Read and become familiar with all instructions, warnings, and cautions before using this device.

A) Site Selection

1. Under direct visualization, examine the patient's cervix to determine the integrity of the cervix.

2. Once the Uterine ElevatOR PRO with OccludOR Balloon™ which has a cervical cup (2, Figure 1) that most closely approximates size of cervix (Figure 2).

3. Ensure the cup size selected is larger than the cervix. This places the vaginal fornices on stretch, allowing the cervix to be contained within the internal diameter of the cup, thus preventing the unintended inclusion of structures such as the uterine in the dissection plane during uterine division and/or colpotomy.

B) Insertion into the Uterine Cavity

1. Remove the Uterine ElevatOR PRO with OccludOR Balloon™ from its sterile packaging and inspect for damage caused by shipping. Discard if damage is noted.

2. Prior to insertion, test the integrity of the intrauterine balloon by drawing 2-3cc of air into a standard fuer lock syringe and insert into the inflation valve (7, Figure 1). Inject air to confirm inflation. Discard the device if the intrauterine balloon does not remain inflated. Discard and use another Uterine ElevatOR PRO with OccludOR Balloon™ unit. After the successful balloon inflation, deflate the balloon by rinsing with the syringe and remove syringe.

3. The patient should be placed in the lithotomy dorsal position.

4. With the cervix under direct visualization, the surgeon should grasp anterior cervix with a simple or similar coated instrument.

5. If necessary, sound the uterus to determine orientation. The gradients are provided as guide for comparison to a standard fuer lock sound.

6. If required, place the cervical cup (2, Figure 1) onto the Uterine ElevatOR PRO with OccludOR Balloon™ shaft.

7. Insert the cervical cup until desired, sutures should be placed prior to the insertion of the Uterine ElevatOR PRO with OccludOR Balloon™.

8. Lubricate the tip and intrauterine balloon with a sterile surgical lubricant.

9. Carefully insert the proximal tip of the Uterine ElevatOR PRO with OccludOR Balloon™ through the cervical os until the device is in the desired position within the uterine cavity. If the uterus is retroverted, rotate the Uterine ElevatOR PRO with OccludOR Balloon™ 180° prior to inserting the intrauterine balloon. Once inserted, rotate the Uterine ElevatOR PRO with OccludOR Balloon™ 180° to align the cervical cup with the vaginal canal.

10. Remove the clamp from the cervix.

11. Using the handle (7, Figure 1) stabilize the external end of Uterine ElevatOR PRO with OccludOR Balloon™ to prevent inadvertent displacement or forward movement of the Uterine ElevatOR PRO with OccludOR Balloon™ until positioned.

12. Reattach the syringe to the pilot balloon (7, Figure 1) and inflate the intrauterine balloon with air until you feel resistance or 10cc of air has been injected, whichever comes first. DO NOT exceed 10cc of air.

C) Positioning the Cervical / Vaginal Cup

1. Move the cervical cup (2, Figure 1) along the Uterine ElevatOR PRO with OccludOR Balloon™ shaft until it seats the outer edge surpasses the vaginal fornix. Ensure proper placement of the cervical cup before proceeding. If lifting the cervical cup, please proceed by gently pulling on the cervical cup.

2. Slide the vaginal cup (3, Figure 1) along the Uterine ElevatOR PRO with OccludOR Balloon™ shaft until it meets the back side of the cervical cup and is completely seated against the vaginal walls to prevent air leak and to prevent inadvertent displacement of the balloon to backward pressure (Figure 4).

3. Move the locking assembly (4, Figure 1) and secure in place by turning the thumbscrew (5, Figure 1) clockwise and slide back towards the handle (6, Figure 1).

4. Apply gentle traction to the Uterine ElevatOR PRO with OccludOR Balloon™ shaft to ensure that device is secured and engaged to the uterus.

D) OccludOR™ Balloon

1. Test inflate the OccludOR™ Balloon (Figure 5 and 6) with 20-60 cc of sterile saline prior to patient insertion in order to verify proper operation of the occluding balloon and then deflate prior to patient insertion.

2. Prior to making colpotomy incisions, inflate the OccludOR™ Balloon with sufficient saline or water to cause the inflated balloon to make continuous radial contact with the vaginal wall. 90 cc of inflation will generally be sufficient.

E) Removal and Disposal of Uterine ElevatOR PRO with OccludOR Balloon™

1. Reinflate an empty syringe to the inflation valve at the end of the pilot balloon (7, Figure 1); fully withdraw air from the intrauterine balloon to deflate.

2. Unlock the locking assembly (4, Figure 1) by turning the thumbscrew (5, Figure 1) anti-clockwise and slide back towards the handle (6, Figure 1).

3. Separate the tissue from the vaginal cup using fingers.

4. Fully release the vaginal cup (3, Figure 1) and slide back (Figure 1).

5. Gently remove the device from the uterus. DO NOT use excessive force to avoid traumatizing the vaginal canal.


6. Before disposing of the Uterine ElevatOR PRO with OccludOR Balloon™, visually inspect the device to check it is intact and all components have been retrieved.

7. Upon completion of the procedure, the Uterine ElevatOR PRO with OccludOR Balloon™ should be discarded in biohazardous waste as per standard procedure.

Intrauterine Injections

The Uterine ElevatOR PRO with OccludOR Balloon™ has a channel for injection of appropriate dyes or marking materials. Remove the cap and insert a luer lock syringe in the dye injection port at the rear of the handle. Use slow and steady pressure for injection of prescribed fluids. Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Figure 1

1. Intrauterine Balloon
2. Cervical Cup
3. Vaginal Cup
4. Locking Assembly

5. Thumbscrew
6. Handle
7. Inflation Valve with Balloon

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