COROLENE®

CE 0459

Français

English

Español

Italiano

Deutsch

Nederlandse

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Polska

Română

Ceská

Magyar

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ENGLISH	INSTRUCTION FOR USE	ET1090 (V14) – 10/2021
	COROLENE® CE 0459 NON ABSORBABLE SURGICAL SUTURE, POLYPROPYLENE MONOFILAMENT	
DESCRIPTION	The synthetic non-absorbable surgical suture COROLENE® are composed of an isotactic polypropylene monofilament, with or without PLEDGET or stop points of different dimensions. COROLENE® suture are white dyed with Titanium dioxide (Ci 77891) or blue dyed with copper Phthalocyanine (Ci 74160) to enhance visibility. The COROLENE® suture comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for non-absorbable surgical sutures. COROLENE® suture elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. The synthetic non-resorbable surgical suture COROLENE® are available in various lengths, diameters (USP/EP) and quantities with surgical needles.	Side effects associated with the use of this medical device include: - wound dehiscence, - calculi formation in urinary or biliary tractus when prolonged contact with salt solutions occurs, - sinus formation - Risk, if existing microbiological infection, that suture, as any foreign substance, acute infection, - Minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body, - Transitory local irritation. Some complications may result from the surgical technical itself.
INDICATIONS	COROLENE® suture are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular and vascular surgery, in ophthalmic surgery, in plastic surgery and in neurological surgery. COROLENE® suture can be used for laparoscopic surgery and abdominal aorta surgery.	VALIDITY / STORAGE CONDITIONS Do not use after the expiry date. Recommended storage conditions: must be stored in the original packaging, below 40°C, in a dry place, away from light.
CONTRAINDICATIONS	Inflammatory and ruptured abdominal aortic aneurysms are contraindications for total laparoscopic repair.	SYMBOLES USED ON LABELLING = PLEDGET = Stop point = Medical device = Unique Device Identifier = Do not reuse = Do not sterilize = Use by date (YYYY-MM-DD) = Date of manufacture (YYYY-MM-DD) = Manufacturer = Batch code = Do not use if package is damaged = Keep away from sunlight = Upper temperature limit = Sterilized using ethylene oxide = Single sterile barrier system = Catalogue number = Caution = Consult instructions for use = CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC
PRECAUTIONS	In order to respect the integrity of the suture material (thread and PLEDGET if need be), care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders (especially with laparoscopic surgery). This device is thermo fused. Pay particular attention while using laser or electro-cautery. During the surgery, grasp the needle 2/3 of the body from the point. If the needle is hold too close to the needle point, penetration qualities should be affected. If it is hold too close to the hole the needle can be twisted and broken at this point. Changing the bending of the needle leads to its breakage. Do not reuse, reprocess, or sterilize Reusing, reprocessing or sterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death. Reusing, reprocessing or sterilizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness. For single patient use only - Contents supplied sterile Do not use if the product or the pouch appears to be damaged.	All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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