SINUSORB® PGA

Français (French)	Suture chirurgicale résorbable, tresse en acide polyglycolique enduite ou monofil en acide polyglycolique non enduit	2
English	Absorbable surgical suture, coated polyglycolic acid braid or uncoated polyglycolic acid monofilament	3
Español (Spanish)	Sutura quirúrgica absorbible, trenza de ácido poliglicólico recubierto o monofilamento de ácido poliglicólico sin recubrimiento	4
Italiano (Italian)	Italano (Italian)	5
Polski (Polish)	Nit chirurgiczna wchłaniająca, plecionka powłokana z kwasu poliglikolowego lub monofilamentowa niepowłokana z kwasu poliglikolowego	6
Deutsch (German)	Resorbierbares chirurgisches Nahtmaterial, beschichteter Geflochtener (Polyglykolsäurefaden) oder unbeschichteter monofiler Polyglykolsäurefaden	7

PETERS SURGICAL
Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT – FRANCE
TEL : +33 (0)1 48 10 62 62 - FAX : +33 (0)1 48 91 22 99
www.peters-surgical.com
E-MAIL: peters@peters-surgical.com

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Marquage CE initial : 2005

SINUSORB® PGA C0459 ABSORBABLE SURGICAL SUTURE, COATED POLYGLYCOLIC ACID Braid OR UNCOATED POLYGLYCOLIC ACID MONOFILAMENT

DESCRIPTION
The sterile synthetic absorbable surgical sutures SINUSORB® PGA are composed of homopolymer of glycolic acid either coated with polyacrylate and calcium stearate or not coated.
The structure of synthetic absorbable surgical sutures SINUSORB® PGA is as follows :
- coated braids for 2 USP 9/0 (0,4 dec).
- uncoated monofilament for USP 9/0 (0,3 dec) and USP 10/0 (0,2 dec).
The sutures SINUSORB® PGA and the coating are inert, noncollagenous, and nonantigenic.
The synthetic absorbable surgical sutures SINUSORB® PGA are available in dyed violet color (D&C violet n°2, CI 60725) or undyed.
The synthetic absorbable surgical sutures SINUSORB® PGA comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for synthetic absorbable surgical sutures.
Expected clinical performances of the synthetic absorbable surgical sutures SINUSORB® PGA are a closure completion without dehiscence and an efficient wound healing without infection.
The synthetic absorbable surgical sutures SINUSORB® PGA elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.
Progressive loss of tensile strength and eventual absorption of the suture occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and metabolized in the body.
The absorption process begins as a loss of tensile strength followed by a loss of mass.
Absorbable sutures SINUSORB® PGA retain in vivo approximately 65 % of their original tensile strength at two weeks post implantation, and more than 40% at three weeks for medium diameter sutures and 30% for small diameter sutures.
The absorption of the sutures SINUSORB® PGA is essentially complete between 60 and 90 days.
The sterile synthetic absorbable surgical sutures SINUSORB® PGA are available in various lengths, diameters (USP/EP) and quantities, single or double armed surgical needles.

INDICATIONS
The synthetic absorbable surgical sutures SINUSORB® PGA are indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic surgery, but not for use in cardiovascular and neurological surgery.

CONTRAINDICATIONS
The use of these sutures is contraindicated in patients with known sensitivities or allergies to their components.
These sutures, being absorbable, should not be used where extended approximation of tissue is required.
Use of absorbable suture may be unsuitable in elderly, malnourished patients, or in patients suffering from cancer, anemia, obesity, diabetes, infection or other conditions which may delay wound healing and especially on sites where expansion, stretching or distension occur.

PRECAUTIONS
Under some circumstances, notably orthopaedic procedure, immobilization by external support may be employed at the discretion of the surgeon.
Skin sutures which must remain in place longer than 7 days may cause localized irritation, they should be removed if necessary.
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
The absorbable sutures are thermofusible. Pay particular attention while using laser or electro-cautery.
During the surgery, grasp the needle 2/3 of the body from the point.
If the needle is held too close to the needle point, penetration qualities should be affected. If it is held too close to the hole the needle can be twisted and broken at this point. Changing the bending of the needle leads to its breakage.
Knot security requires the standard surgical techniques of square knots according to surgical circumstances and the surgeon's experience. Use of additional loops may be particularly appropriate when knotting monofilaments.
In case of oversewing, the surgeon will use non slipknots, self-locking knot at both ends.
Do not use if the product or the pouch appears to be damaged.
Apply a strong pull on the suture to verify the thread strength.
For single patient use only - Contents supplied sterile
Do not reuse, reprocess, or restertize. Reusing, reprocessing or restertizing may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reusing, reprocessing or restertizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious diseases) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
The technique, method and products used depend on clinical experience and are left to the surgeon's discretion in order to ensure the best clinical results.
Users should be familiar with surgical procedures and techniques involving synthetic absorbable sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.
The use of supplemental non absorbable sutures should be considered by surgeon in the closure of sites which may undergo expansion, stretching, or distension, or which may require additional support.
Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness.

USE
Inspect the packaging to confirm its integrity before opening.
Open the peelable pouch and take the tearable aluminium pouch.
Tear the pouch by gripping the end at the right side of the arrow and draw it towards oneself.
Access to the needle is direct, grab it and pull the suture out.

ADVERSE EFFECTS
Adverse effects associated with the use of this device include:
- wound dehiscence,
- wound supuration and bleeding with sinus formation,
- risk, if existing microbiological infection, that suture, as any foreign body, accentuate infection,
- moderate tissue inflammatory response, characteristic of foreign body response,
- localized irritation when skin sutures are left in place for greater than 7 days,
- calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs,
- transitory local irritation,
- eversion,
- scars stretching.
Some complications may result from the surgical technical itself.

VALIDITY / STORAGE CONDITIONS
Do not use after the expiry date.
The product must be stored in its original packaging, at a temperature lower than 25°C, in a dry place, and in the light.

SYMBOLES USED ON PACKAGING

	= Medical device
	= Unique Device Identifier
	= Do not reuse
	= Do not restertize
	= Use by date (YYYY-MM-DD)
	= Date of manufacture (YYYY-MM-DD)
	= Manufacturer
	= Batch code
	= Do not use if package is opened or damaged
	= Keep dry
	= Keep away from sunlight
	= Upper limit of temperature
	= Sterilized using ethylene oxide
	= Single sterile barrier system with protective packaging inside
	= catalogue number
	= Caution
	= Consult instructions for use
	= (for USA) CAUTION: US Federal Law restricts this device to prescription only
	= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

PETERS SURGICAL
Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT - FRANCE
TEL : +33 (0)1 48 10 62 62 - FAX : +33 (0)1 48 91 22 99
www.peters-surgical.com
E-MAIL: peters@peters-surgical.com

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