POLYPROPYLEN, SOIE/SILK, NYLON

SUTURES NON RESORBABLES, CHIRURGIE OPHTHALMIQUE

DESCRIPTION

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Marquage CE initial NYLON : 1997

POLYPROPYLEN, SOIE/SILK, NYLON

NON ABSORBABLES SUTURES, OPHTHALMIC SURGERY

PRODUCT

Sterile needle/suture combinations, double packed in two pouches consisting of:

- non absorbable monofilament or multifilament suture,
- needle made of stainless medical steel, each single or double armed.

The following suture materials are being attached:

Polypropylen	Extruded monofilament of polypropylen
	Blue color
Soie / Silk	Made of spun of the silk worm, waxed multifilament
	Black color
Nylon	Extruded monofilament of polyamide 6.6
	Black color

INDICATIONS

POLYPROPYLEN, SOIE/SILK, NYLON sutures are intended for use in ophthalmic surgery.

ACTIONS

POLYPROPYLEN, SOIE/SILK, NYLON sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

As with any polyamide and silk sutures, though non absorbable, gradual loss of tensile strength can occur over a long time in vivo.

CONTRAINDICATIONS

All suture materials can cause a minimal inflammatory reaction.

Tensile strength can be lost over time in the case of nylon due to a progressive hydrolysis and in the case of silk due to degradation in vivo. They should not be used where permanent retention of tensile strength is needed and should only be used where according to surgical practice a non-absorbable suture material is being used.

The use of silk suture is contraindicated in patients with known sensitivities or allergies to silk.

Do not use in direct contact with the heart, the central circulatory system, or the central nervous system.

PRECAUTIONS

As with any suture, care should be taken in handling to avoid damage. Avoid thread crushing or crimping damage due to application of surgical instrument, such as needle holders or forceps.

Beware when laser or electro-cautery is used as the sutures can be thermo fused.

Sutures should be handled preferably with a needle holder at a third or half of the distance between the connection to the braid and the point to avoid penetration difficulties, rupture or needle bending. Changing the bending of the needle leads to a torsion resistance loss and can cause the breakage.

Knot security requires the standard surgical techniques of square ties according to surgical circumstances and the surgeon's experience. Use of additional throws may be particularly appropriate when knotting monofilaments.

The design of this device doesn't allow its reuse.

Do not sterilize. Sterile unless packaging has been opened or damaged. Destroy after use.

Users should be familiar with surgical procedures and techniques using non-absorbable sutures employing **POLYPROPYLEN, SOIE/SILK, NYLON**, as the risk of wound dehiscence may vary with the site of implantation and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness.

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

Products should be destroyed in appropriate incinerator in respect to local and national hospital wastes regulation.

SIDE EFFECTS

Adverse effects associated with the use of this medical device include:

- wound dehiscence;
- calculi formation in urinary or biliary tractus when prolonged contact with salt solutions occurs;

- risk, if existing microbiological infection, that suture, as any foreign substance, accentuate infection;
- minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body, and transitory local irritation.

Allergic response in patients that are known to be sensitive to silk has been observed.

SUPPLY FORMS

The non-absorbable surgical sutures **POLYPROPYLEN, SOIE/SILK, NYLON** are available in various lengths, diameters (USP/EP) and quantities. Each box contains 12 single or double armed sutures.

SHELF LIFE / PREVENTION

Recommended storage conditions: below 40°C, away from moisture and direct heat. Do not use after expiration date.

SYMBOLS USED ON LABELLING

= Medical device

= Do not reuse

= Do not resterilize

= Use by date (YYYY-MM-DD)

= Date of manufacture (YYYY-MM-DD)

= Manufacturer

= Batch number

= Do not use if package is opened or damaged

= Keep away from sunlight

= Keep dry

= Upper temperature limit

= Sterilized using ethylene oxide

= Single sterile barrier system

= Single sterile barrier system with protective packaging inside

= catalogue number

= Unique Device Identifier

= Caution

= Consult instructions for use

= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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