**GLACTIME® CE0459****SUTURE CHIRURGICALE RESORBABLE, TRESSE DE PGLA ENDUITE****DESCRIPTION**

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GLACTIME® CE0459**ABSORBABLE SURGICAL SUTURE, COATED BRAIDED PGLA****DESCRIPTION**

The sterile synthetic absorbable surgical sutures GLACTIME® are composed of copolymer of glycolic and lactic acids and coated with copolymer of glycolic and lactic acid and calcium stearate.

The sutures and coating are inert, nonreactive and nonantigenic.

The synthetic absorbable sutures GLACTIME® are available in dyed violet color (D&C violet n°2, CI 60725) or undyed.

The synthetic absorbable sutures GLACTIME® comply with the requirements of the European Pharmacopeia (EP) for synthetic braided absorbable surgical sutures and the United States Pharmacopeia (USP), except for minor size in diameter.

Expected clinical performance of the synthetic absorbable surgical sutures GLACTIME® is to provide absorption without desiccation and an efficient dissolution process begins with a loss of tensile strength followed by a loss of mass.

Absorbable sutures GLACTIME® retain in vivo approximately 65 % of their original strength at week 3 post implantation, and more than 40 % at three weeks.

Absorption of the sutures GLACTIME® is essentially complete between 56 and 70 days.

The synthetic absorbable sutures GLACTIME® are available in various lengths, diameters (USP/EP) and quantities with or without surgical needles.

In general, the synthetic absorbable surgical sutures GLACTIME® are indicated for use in general tissue approximation and/or ligation, including use in ophthalmic surgery.

CONTRADICTIONS

Absorbable sutures GLACTIME® must not be used in cardiovascular and neurological tissues.

The use of these sutures is contraindicated in patients with known sensitivities or allergies to their components.

These sutures, being absorbable, should not be used where extended periods of time are required.

Use of absorbable sutures GLACTIME® is not suitable in elderly or malnourished patients, or in patients suffering from cancer, anaemia, obesity, diabetes, infection or other conditions which may delay wound healing and especially on sites with expansion, stretching or delamination.

PRÉCAUTIONS

Under some circumstances, notably orthopaedic procedure, immobilization or external support may be employed at the discretion of the surgeon.

Suturing which remains too long can then 7 days may cause localized irritation, they should be removed if necessary.

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to use of sharp instruments.

The absorbable sutures GLACTIME® are thermoresistant. Pay particular attention while using laser or electro-cautery.

If the needle is held too close to the needle point, penetration qualities should be affected, it should be held close to the hole the needle can become localized irritation, they should be removed if necessary.

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to use of sharp instruments.

Knot security requires the standard surgical techniques of square knots according to surgical experience and the surgeon's experience. Use of additional loops may be particularly appropriate.

In case of over-tensioning, the surgeon may use a self-locking siphon knot, according to surgical experience and the surgeon's experience.

Do not use if the product has been punctured or damaged.

For single patient use only - Contents supplied sterile.

No reprocessing, resterilize, reprocessing, or resterilizing may compromise the structural integrity of the device and/or cause a risk of contamination of the device or cause patient infection or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device or cause patient infection or death. Reusing, reprocessing or resterilizing may also create a risk of transmission of an infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

The technique, method and products used depend on clinical experience and left to the surgeon's discretion in order to ensure the best clinical results.

Users should be familiar with surgical procedures and techniques involving synthetic absorbable sutures for wound closure, as risk of wound

VALIDITY / STORAGE CONDITIONS

Do not use after the expiry date.

Do not use in the original packaging, at a temperature lower than 25°C in a place of light.

SYMBOLS USED ON PACKAGING

= Medical device

= Unique Device Identifier

= Do not reuse

= Do not sterilize

= Use by date (YYYY-MM-DD)

= Date of manufacture (YYYY-MM-DD)

= Manufacturer

= Batch code

= Do not use if package is opened or damaged

= Keep dry

= Keep away from sunlight

= Upper limit of temperature

= Sterilized using ethylene oxide

= Single sterile barrier system with protective packaging inside

= catalogue number

= Caution

= Consult instructions for use

= CE Mark. Product conform to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

PRECAUTIONS

These sutures are accidently caused with the site of application and the suture material

desiccance may follow with a loss of tensile strength followed by a loss of mass.

The sutures GLACTIME® are composed with respect to drainage and closure or contortion of the suture or which may necessitate a support additional.

Use of sutures GLACTIME® and the mélange d'enduction are inert, privé de collagène ou non antigenic.

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