CARDIOFLON®

Evolution

CE 0459

Français	Suture chirurgicale non résorbable, tresse de polyester enduite
English	Non absorbable surgical suture, coated polyester braid
Español	Sutura quirúrgica no absorbible, hilo recubierto de poliéster
Italiano	Sutura chirurgica non assorbibile, treccia di poliestere rivestito
Deutsch	Nicht resorbierbares chirurgisches Nahtmaterial, beschichteter, geflochtener Polyesterfaden
Nederlandse	Niet-resorbeerbare hechtdraad voor chirurgisch gebruik, strengen gemaakt van gecoat polyester
Ελληνικο	Μη απορροφησιμο χειρουργικό ραμμα, πλέγμα επικαλυμμένο με πολυεστέρα
Polska	Niewchłanianalne nici chirurgiczne z powlekanego splotu poliestrowego
Română	Fir chirurgical de sutură neabsorbabil, din fibre de poliester
Ceská	Nevstřebatelná chirurgická nit, pletený potahovaný polyester
Magyar	Nem felszívódó sebészeti varrófonal, bevonatos poliészter fonál
Русский	Нерассасывающийся хирургический шовный материал, плетеная полиэфирная нить с покрытием
عربي	خيطة جراحي غير قابل للامتصاص، جديلة بوليستر مشمعة
中文	非吸收性外科缝线，聚酯涂层织带（商品名：卡迪隆升级）

PETERS SURGICAL

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Marquage CE initial : 2013

CARDIOFLON® Evolution CE0459

NON ABSORBABLE SURGICAL SUTURE, COATED POLYESTER BRAID

DESCRIPTION

The synthetic non-absorbable surgical sutures **CARDIOFLON® Evolution** are composed of polyester [poly(ethylene terephthalate)] braids, coated, with or without **PLEDGES** of different sizes.

CARDIOFLON® Evolution sutures are green dyed (D&C green n°6 CI 61565) and white (undyed). Some packs are composed of green and white braids to facilitate location during use.

CARDIOFLON® Evolution sutures comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for non-absorbable surgical sutures, except for some diameters.

The synthetic non-absorbable surgical sutures **CARDIOFLON® Evolution** are available in various lengths, diameters (USP/EP) and quantities with or without needles.

CARDIOFLON® Evolution sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

INDICATIONS

CARDIOFLON® Evolution sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, vascular and ophthalmic surgeries.

CONTRAINDICATIONS

None known.

PRECAUTIONS

In order to respect the integrity of the suture material (thread and **PLEDGET** if need be), care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to application of surgical instruments such as forceps or needle holders (especially with laparoscopic surgery).

This device is thermo fused. Pay particular attention while using laser or electro-cautery.

During the surgery, grasp the needle 2/3 of the body from the point. If the needle is held too close to the needle point, penetration qualities should be affected, if it is held too close to the hole the needle can be twisted and broken at this point. Changing the bending of the needle leads to its breakage.

Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death.

Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness.

For single patient use only - Contents supplied sterile

Do not use if the product or the pouch appears to be damaged.

USE

Inspect the packaging to confirm its integrity before opening.

As any suture material, the technique (knots safety), methods and associated products must be adapted by the surgeon according to the clinical context and his experience to ensure patient's safety. Additional knots can be added for monofilament.

Users should be familiar with surgical procedures and techniques involving synthetic non-absorbable sutures, as risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

This product can be used only by or on the order of a physician.

Any device which has been used or partially used must be destroyed in proper incinerator in accordance with the local or national regulations concerning the disposal of hospital waste.

SIDE EFFECTS

Side effects associated with the use of this medical device include:

- wound dehiscence,
- calculi formation in urinary or biliary tractus when prolonged contact with salt solutions occurs,
- Sinus formation
- Risk, if existing microbiological infection, that suture, as any foreign substance, accentuate infection,
- Minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body,
- Transitory local irritation.

Some complications may result from the surgical technical itself.

VALIDITY / STORAGE CONDITIONS

Do not use after the expiry date.

Recommended storage conditions: must be stored in the original packaging, below 40°C, in a dry place, away from light.

SYMBOLS USED ON LABELLING

	= PLEDGET
	= Medical device
	= Unique Device Identifier
	= Do not reuse
	= Do not resterilize
	= Use by date (YYYY-MM-DD)
	= Date of manufacture (YYYY-MM-DD)
	= Manufacturer
	= Batch code
	= Do not use if package is damaged
	= Keep away from sunlight
	= Keep dry
	= Upper temperature limit
	= Sterilized using ethylene oxide
	= Single sterile barrier system
	= Catalogue number
	= Caution
	= Consult instructions for use
	= (for USA): CAUTION: US Federal Law restricts this device to prescription only.
	= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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