IFABOND®

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EN

SYNTHETIC SURGICAL GLUE

Before using the product, check the entire package (do not use if the blister is damaged). Read all the instructions.

SINGLE USE PRODUCT - DO NOT RE-USE - DO NOT RE-STERILISE - DO NOT USE IF THERE IS EVIDENCE OF DAMAGE**DESCRIPTION**

IFABOND® surgical glue is a class III absorbable and sterile medical device. The glue is a translucent fluid, made from n-hexyl cyanoacrylate monomers. IFABOND® is packaged in a hermetically sealed polypropylene vial. Each vial contains 0.5 ml, 1 ml or 1.5 ml of glue and is packed in a double sterile blister, supplied with a syringe and a sterile, single-use intra-muscular needle for extracting the glue from the vial.

INTENDED USE

IFABOND® is intended to be used during surgical procedures, in open surgery and in laparoscopic surgery, for its adhesive and hemostatic sealant actions.

INDICATIONS

IFABOND® surgical glue is intended to be used during open surgery and laparoscopic surgery.

These surgeries for the adult population include digestive and visceral surgery for the treatment of hernia, bariatric surgery for sleeve-bypass, urological and gynecological surgery with sacrocolpopexy and breast surgery for mastectomy and axillary curage. The glue can be used in the pediatric population for circumcision.

TARGET POPULATION

IFABOND® product can only be used by or on the order of a physician. Users should be familiar with surgical procedures and techniques involving synthetic glue. Globally, the use of glue requires a specific training provided to surgeon and operating room nurse. The device is intended for an adult population and for pediatric population for circumcision.

ACTION MODE

IFABOND® mode of action is based on adhesion and hemostatic sealant action promoting reaction in tissues, which is followed by gradual encapsulation of the glue by fibrous connective tissue and progressive medical device degradation.

IFABOND® may be applied alone or in combination with sutures. This medical device is designed to offer an alternative to mechanical stapling, sutures and any other means used during surgical procedures. It therefore makes it possible not only to close superficial or deeper wounds without tension on them but also to fit a prosthesis to muscles, tendons or other organs *in vivo*.

CLINICAL PERFORMANCE AND BENEFITS

The Clinical performances associated to IFABOND® are non-traumatic fixation (no possible transfixion) for tissue or mesh, tissues healing and closure of superficial or deeper wounds without tension.

The Patients' benefits demonstrated after clinical evaluation are shorter operating time, low length of hospital stay, less postoperative pain, shorter time to return to daily activities, low risk of recurrence.

CLINICAL SAFETY

The safety of patients and users is guarantee when IFABOND® is used. To date, no side effects related to the device are known, except for a risk related to the use of glue in significant quantity producing an exothermic reaction. The local inflammation caused by this reaction does not call into question the safety and health of the patient.

CONTRA-INDICATIONS

- The glue should not be applied in direct contact with fragile and/or damaged tissue which might not withstand thermal, inflammatory or mechanical stresses.
- The glue should not be applied in direct contact with cerebral tissue.
- Use of the glue is contra-indicated in ophthalmic and orthopaedic surgery.
- The glue is not suitable for intra-vascular use.
- The glue should not be applied to wounds, incisions or openings subject to high levels of tension, unless in conjunction with other securing and closing methods.
- Use of the glue is contra-indicated for any infected wound or where there is a risk of gangrene.
- The glue is not intended for use in a cutaneous application
- The glue should not be used on pregnant women, on patients with a known hyper-sensitivity to cyanoacrylate or formaldehyde or on patients with a systemic pre-operative infection, with uncontrolled diabetes, or with conditions known to interfere with the healing process.

PRECAUTIONS

Apply the glue in fine layers or drop by drop (approximately ratio one drop per square centimeter), to obtain maximum adhesion and good scarring. It is important not to apply several drops of glue to the same point. Do not use if the glue cannot be applied in a controlled manner (quantity, flow, and system of delivery). Unless different medical prescription, a small amount of glue should be applied respecting the ratio of 1 drop per cm².

If too great a quantity is accidentally applied, remove it immediately with a damp pad/swab or rinse with saline solution

Indeed, an excess of polymerized product can form rigid accumulations of glue which should be removed to prevent any risk of trauma to tissue or organs in the area around the application zone.

WARNINGS

IFABOND® should not be used together with devices or instruments made of polycarbonate.

IFABOND® is ready for use and should not be diluted or mixed with other substances.

IFABOND® is a rapid polymerization adhesive capable of joining body tissue to many other materials such as latex or steel. Accidental contact with any other tissue, surfaces or equipment should be prevented and avoided, if they cannot be copiously rinsed with lukewarm water or easily cleaned using a solvent such as acetone. In the event of accidental contact with the eyes or skin, rinse abundantly with water.

If applied unintentionally, do not attempt to unstick IFABOND®, unless it is absolutely necessary, as this could cause damage to the tissue. The glue should always be applied with perfect control over the amount laid down. Excess glue could cause a slight release of heat and poor adhesion.

Check that all items used in conjunction with the glue are sterile and compatible, to avoid any deterioration of its properties.

Using IFABOND® in cardiovascular surgery, maxillofacial and neurosurgery are currently not supported by sufficient clinical data.

APPLICATION

- Open the extern blister.
- Remove the sterile inner blister.
- Place the inner blister directly on the operating table in a sterile environment.
- Open the inner blister and remove the sterile vial.
- Check that the glue is fluid and transparent. The product should be used without any additional preparation.
- After application, the rest of the product should be disposed of with its packaging and should not be used in any other procedure.

During application, the glue remains liquid and then, upon contact with living tissue or damp and/or alkaline mediums, it quickly polymerizes and provides a good mechanical bond. Once polymerization is complete, the glue no longer has any adhesive power. Under normal conditions of use, the glue is resorbed by hydrolytic degradation, the resorption time varying according to the type of tissue and the quantity of glue applied.

Applicator use

- Unscrew the IFABOND® vial's cap and extract the product by using the sterile Luer lock syringe and intra-muscular needle included. Do not purge needle and syringe after collection.
- Remove the needle used to puncture the vial and fit the recommended application accessories or catheter on the syringe without expelling the air in the syringe
- Your IFABOND® applicer device is ready to use.
- For maximum adhesion and safety, the glue should be applied in small quantities. The glue can be applied as a uniform line or drop to drop (≈ 1 drop / cm^2) according to the applicator used.
- Given the low viscosity of the glue (viscosity close to that of water), it is important to apply the right amount of pressure on the syringe piston, in case of drop-by-drop dispensing: excess polymerized product can cause the formation of rigid accumulations of glue which must then be removed.
- It is not necessary to apply several drops of product to the same spot, as this could cause the formation of a more rigid accumulation of glue or result in a longer polymerization time and reduce the quality of the adhesion (but if multiple drops are applied, ensure that the initial ones polymerize completely before new application on the same area).
- As a reference, 1 ml of IFABOND® covers a surface of approximately 50 cm^2 .

Spray use

- Unscrew the IFABOND® vial's cap and insert the vial into an IFAJET® sprayer, ref. VA01.
- Follow the instructions given in the manual supplied with that device IFAJET®.
- In some cases, it is advisable to protect the areas around the area to be sprayed.
- If it proves necessary to apply several layers, wait until the previous layer has completely polymerised; however, this will reduce the flexibility of the film formed and could cause a rigid mass of glue to form.
- As a guide, 1.5 ml of sprayed IFABOND® covers a surface of approximately 150 cm^2 .

STORAGE

IFABOND® must be stored in its original packaging away from direct sunlight and humidity at a temperature of 5°C or below.

VALIDITY

Under recommended storage conditions, IFABOND® has a shelf life of two years (date indicated on the packaging).

Do not use the product after the expiry date.

PRODUCT RANGE

Boxes	Name	References
6 sterile blister packs	0.5 ml vial of IFABOND® glue	IB05
6 sterile blister packs	1 ml vial of IFABOND® glue	IB
6 sterile blister packs	1.5 ml vial of IFABOND® glue	IB+

SINGLE USE DEVICE

IFABOND® is a single-use sterile medical product for use on a single patient only. If re-used, the risks involved are those relating to a breakdown of sterility and cross-contamination.

The product should not be re-sterilised and should not be used if it is damaged. Check the integrity of the package prior to use.

STERILISATION METHOD

The glue, syringe and needle are sterilised in the double blister by beta radiation.

SYMBOLS USED

	Medical device		Batch code		Upper temperature limit
	Catalogue number		Use-by date (YYYY-MM-DD)		Keep dry
	Unique Device Identifier		Do not use if the packaging is damaged		Store away from sunlight
	Date of manufacture (YYYY-MM-DD)		Do not reuse		Caution
	Manufacturer		Do not resterilize		Consult instructions for use
	Sterilized using irradiation		Double sterile barrier system		Product complies with the requirements of directive 93/42/EEC

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