Purefill®

| | | |
|----------------------------|--|----|
| Français (French) | AGENT D'EMBOLISATION LIQUIDE | 2 |
| English | LIQUID EMBOLIC AGENT | 3 |
| Español (Spanish) | AGENTE DE EMBOLIZACIÓN LÍQUIDO | 4 |
| Deutsch (German) | FLÜSSIGES EMBOLISATIONSMITTEL | 5 |
| Português (Portuguese) ... | AGENTE DE EMBOLIZAÇÃO LÍQUIDO | 6 |
| Italiano (Italian) | AGENTE DI EMBOLIZZAZIONE LIQUIDO | 7 |
| Ελληνικά (Greek) | ΥΓΡΟΣ ΠΑΡΑΓΟΝΤΑΣ ΕΜΒΟΛΙΣΜΟΥ | 8 |
| Svenska (Swedish) | FLYTANDE EMBOLIMEDEL | 9 |
| Dansk (Danish) | FLYDENDE EMBOLISK MIDDEL | 10 |
| Norks (Norwegian) | FLYTENDE EMBOLISK MIDDEL | 11 |
| Nederlands (Dutch) | VLOEIBAAR EMBOLISATIEPRODUCT | 12 |
| Русский (Russian) | ЖИДКОЕ СРЕДСТВО ДЛЯ ЭМБОЛИЗАЦИИ | 13 |
| عربية (Arabic) | عامل الانصمام السائل | 14 |
| Polski (Polish) | CIEKŁY ŚRODEK DO EMBOLIZACJI | 15 |
| Magyar (Hungarian) | EMBÓLIA ELLENI FOLYÉKONY SZER | 16 |
| Türkçe (Turkish) | SIVI EMBOLİ ÜRÜNÜ | 17 |
| 中文 (Chinese) | 液体栓塞剂 | 18 |
| Romană (Romanian) | AGENT DE EMBOLIZARE LICHID | 19 |
| Ceská (Czech) | TEKUTÝ EMBOLIZAČNÍ PROSTŘEDEK | 20 |

PETERS SURGICAL

Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT - FRANCE
www.peters-surgical.com - E-MAIL: peters@peters-surgical.com
TÉL.: +33 (0)1 48 10 62 62 - FAX: +33 (0)1 48 91 22 99

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CE 0459

DESCRIPTION

PUREFILL® is a sterile, translucent, liquid embolic agent composed of n-hexyl-cyanoacrylate monomers.

PUREFILL® is packaged in a hermetically sealed polypropylene vial. The vial containing 1 ml of product is packaged in a sterile blister pack. When applied, the embolic agent remains liquid until it comes in contact with tissue or moist environments and then it polymerises in a few seconds. Once polymerisation is completed the embolic agent has no more adhesive power.

PUREFILL® is shipped in an isotherm box containing 2 boxes of 6 blisters.

INDICATIONS

PUREFILL® is a liquid embolic agent indicated in the embolisation of arteriovenous malformations (AVM).

PUREFILL® is a class III implantable medical device which is described as an arterial and venous embolisation implant due to its haemostatic action, and is applicable in the fields of neuroradiology and interventional radiology. This medical device is applied in solution with a contrast liquid.

CONTRAINDICATIONS

Do not use **PUREFILL®** on patients with a history of known hypersensitivity to cyanoacrylate or formaldehyde.

The use of **PUREFILL®** is contraindicated when one of the following conditions exists:

- When optimum positioning of the catheter is not possible
- When occlusion tests indicate an intolerance to the occlusion procedure
- When a vasospasm stops blood flow.
- When significant blood flow prevents the injection of an embolic agent.

PRECAUTIONS

Embolisation to block blood vessels is a very high-risk endovascular procedure. **PUREFILL®** should be used by physicians with training in interventional neuroradiology or peripheral interventional radiology and in-depth knowledge of the disease to be treated, and of angiographic and embolisation techniques.

Ensure that the catheters and ancillary articles used are clean and compatible with the medical device and do not trigger any premature polymerisation (please refer to the instructions for use)

Like all cyanoacrylates, **PUREFILL®** should be used with appropriate syringes (PE or PP): do not use polycarbonate syringes.

Do not use if the device or the blister appears to be damaged.

For single use only - Contents supplied sterile

Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/ or lead to device failure, which in turn may result in patient injury, illness or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PUREFILL® is a translucent liquid which should not be used if it has thickened.

The presence of any impurity in the product could result in rapid polymerisation preventing its reuse.

This is a rapid polymerisation cyanoacrylate. Therefore, prevent and avoid accidental contact with all tissue, surfaces or equipment that cannot be rinsed abundantly with warm water or easily cleaned with a solvent such as acetone.

Any device which has been used or partially used must be destroyed in accordance with the local or national regulations concerning the disposal of hospital waste.

USE

- Inspect the packaging to confirm its integrity before opening
- Remove the blister pack from its storage place a few minutes before use
- Verify the fluidity and transparency of the product
- Prepare the **PUREFILL®** / contrast fluid mixture in an adapted zone: a sterile glass beaker is recommended
- For all occlusion procedures, mix **PUREFILL®** with iodized oil (such as Lipiodolultrafluid®) in 01:02 proportions
- Homogenize, then withdraw the mixture with a Luer lock syringe for injection of 1 ml to 3 ml
- Before injecting, verify the radio-opacity of the mixture
- Position the micro-catheter according to the arterial embolisation radiological procedure after rinsing with a 5% glucose solution
- Connect the syringe to the micro-catheter and inject the mixture manually under radiographic monitoring
- After injection, rapidly remove the micro-catheter to prevent adhesion of the tip

- After each injection, discard the micro-catheter, all the accessories and the remaining mixture

ADVERSE EFFECTS

The use of n-hexyl-cyanoacrylate can cause a mild exothermic reaction (+2°C) that can result in local inflammation if the quantity of product applied at one point is significant.

VALIDITY / STORAGE CONDITIONS

Do not use after use-by date.

PUREFILL® should be kept in its original packaging protected from direct sunlight and humidity, at a storage temperature below 5°C.

SYMBOLS USED ON PACKAGING

= Medical device

= Manufacturer

= Catalogue reference

= Batch number

= Date of manufacture (YYYY-MM-DD)

= Use until (YYYY-MM-DD)

= Sterilized by Irradiation.

= Single sterile barrier system with protective packaging inside

= Do not reuse

= Do not resterilize

= Do not use if package is opened or damaged

= Unique Device Identifier

= Caution

= Consult instructions for use

= Protect from humidity

= Protect from direct sunlight

= Upper temperature limit

= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC.

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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