ET995 (V12) – 10/2021

CARDIONYL®

CE 0459

Français
English
Español
Italiano
Deutsch
Nederlandse
Ελληνικό
Polski
Română
Česká
Magyar
Русский
عربي
中文

Suture chirurgicale non résorbable, monofil de polyamide
Non absorbable surgical suture, polyamide monofilament
Sutura quirúrgica no absorbible de monofilamento de poliamida
Sutura chirurgicala non absorbibilă în monofilament di poliamidă
Niet-absorbeerbare chirurgische hecht draad, polyamide monofilament
Μη απορροφώσιμο χειρουργικό ραμμά, από μονοκλώριο πολυαμιδιο
Monofilamentowa nici chirurgiczna niewchłaniałana z poliamidu
Fir chirurgical de sutură neabsorbabil, din poliamidă, monofilament
Nevstřebatelný monofilamentní chirurgický šicí materiál z polyamidu
Nem felszívódó, poliamid alapanyagú, monofil sebészeti fonal
Нерассасывающийся монофиламентный хирургический шовный материал на основе полиамида покрытием
خطوط الغرز الجراحية الغير قابلة للامتصاص، مكونة من خياف من البوليستر المزوج بالتقرون
不可吸收外科手术缝线，聚酰胺单丝（商品名：卡皮迪尔）

PETERS SURGICAL
Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT - FRANCE
TEL : +33 (0)1 48 10 62 62 - FAX : +33 (0)1 48 91 22 99
www.peters-surgical.com
E-MAIL: peters@peters-surgical.com

Immeuble AURELIUM - 1 Cours de l'Île Seguin - 92100 BOULOGNE-BILLANCOURT - FRANCE
TEL : +33 (0)1 48 10 62 62 - FAX : +33 (0)1 48 91 22 99
www.peters-surgical.com
E-MAIL: peters@peters-surgical.com

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ENGLISH	INSTRUCTION FOR USE	ET995 (V12) – 10/2021
	CARDIONYL® CE 0459 NON ABSORBABLE SURGICAL SUTURE, POLYAMIDE MONOFILAMENT	
DESCRIPTION	The synthetic non-absorbable surgical suture CARDIONYL® are composed of a polyamide monofilament, with or without PLEDGET of different sizes. CARDIONYL® suture are blue dyed (FD&C blue n°2 CFR 741102). CARDIONYL® sutures comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for non-absorbable surgical suture. CARDIONYL® sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. As with any polyamide suture, though non absorbable, gradual loss of tensile strength can occur over a long time in vivo. The synthetic, non-absorbable surgical suture CARDIONYL® are available in various lengths, diameters (USP/EP) and quantities with surgical needles.	- Sinus formation - Risk, if existing microbiological infection, that suture, as any foreign substance, accentuates infection. Minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body. - Transitory local irritation. Some complications may result from the surgical technical itself. VALIDITY / STORAGE CONDITIONS Do not use after the expiry date. Recommended storage conditions: must be stored in the original packaging, below 40°C, in a dry place, away from light. SYMBOLES USED ON LABELLING = PLEDGET = Medical device = Unique Device Identifier = Do not reuse = Do not resterilize = Use by date (YYYY-MM-DD) = Date of manufacture (YYYY-MM-DD) = Manufacturer = Batch code = Do not use if package is damaged = Keep away from sunlight = Keep dry = Upper temperature limit = Sterilized using ethylene oxide = Single sterile barrier system = Catalogue number = Caution = Consult instructions for use = (for USA): CAUTION: US Federal Law restricts this device to use by prescription only. = CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.
INDICATIONS	CARDIONYL® suture are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, vascular, ophthalmic and neurological surgery.	
CONTRAINDICATIONS	CARDIONYL® suture should not be used with synthetic material prosthesis. Using these suture in pre-term infant population is not supported by sufficient clinical data.	
PRECAUTIONS	In order to respect the integrity of the suture material (thread and PLEDGET if need be), care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to application of surgical instruments such as forceps or needle holders (especially with laparoscopic surgery). This device is thermo fused. Pay particular attention while using laser or electro-cautery. During the surgery, grasp the needle 2/3 of the body from the point. If the needle is hold too close to the needle point, penetration qualities should be affected, if it is held too close to the hole the needle can be twisted and broken at this point. Changing the bending of the needle leads to its breakage. Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and or lead to device failure, which in turn may result in patient injury, illness or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness. For single patient use only - Contents supplied sterile Do not use if the product or the pouch appears to be damaged.	
USE	Inspect the packaging to confirm its integrity before opening. As any suture material, the technique (knots safety), methods and associated products must be adapted by the surgeon according to the clinical context and his experience to ensure patient's safety. Additional knots should be considered with monofilament thread. Users should be familiar with surgical procedures and techniques involving synthetic non absorbable sutures, as risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds. This product can be used only by or on the order of a physician. Any device which has been used or partially used must be destroyed in proper incinerator in accordance with the local or national regulations concerning the disposal of hospital waste.	
ADVERSE EVENTS	Side effects associated with the use of this medical device include: - wound dehiscence, - calculi formation in urinary or biliary tractus when prolonged contact with salt solutions occurs,	

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