PREMIO®

CE 0459

Français	Suture chirurgicale non résorbable, monofil de PVDF
English	Non absorbable surgical suture, PVDF monofilament
Español	Sutura quirúrgica no absorbible, monofilamento de fluoruro de polivilideno
Italiano	Sutura chirurgica non assorbibile in monofilamento di PVDF
Deutsch	Nicht resorbierbares chirurgisches Nahtmaterial, PVDF Monofil
Nederlandse	Niet-absorbeerbare chirurgische hechtdraad, PVDF monofilament
Ελληνικό	Μη απορροφησιμο χειρουργικό ραμμα μονοκλωνο φθοριουχο πολυβινυλιδιο
Polska	Monofilamentowa nić chirurgiczna niewchłaniałna z PVDF
Română	Fir chirurgical de sutură neabsorbabil, PVDF monofilament
Ceská	Monofilamentní nevstřebatelný chirurgický šicí materiál z PVDF
Magyar	Nem felszívódó, PVDF alapanyagú, monofil sebészeti fonal
Русский	Рассасывающийся хирургический шовный материал из монофиламентного поливинилиденфторида
عربي	خيوط الغرز الجراحية غير القابلة للامتصاص، خيط أحادي الألياف من البولي فينيليدين فلورايد
中文	不可吸收外科手术缝线，PVDF 单丝（商品名：派迷欧）

PETERS SURGICAL

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The synthetic non-absorbable surgical sutures **PREMIO®** are composed of a PVDF (polyvinylidene fluoride) monofilament, with or without **PLEDGET** of different sizes.

PREMIO® sutures are blue-violet dyed (CI Solvent blue 104).

PREMIO® sutures comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for non-absorbable surgical sutures, except for some diameters.

PREMIO® sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

The synthetic non-absorbable surgical sutures **PREMIO®** are available in various lengths, diameters (USP/EP) and quantities with surgical needles.

INDICATIONS

PREMIO® sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, vascular and neurological surgery.

CONTRAINDICATIONS

None known.

PRECAUTIONS

In order to respect the integrity of the suture material (thread and **PLEDGET** if need be), care should be taken to avoid damage from handling. Avoid thread crushing or crimping damage due to application of surgical instruments such as forceps or needle holders (especially with laparoscopic surgery).

This device is thermo fused. Pay particular attention while using laser or electro-cautering.

During the surgery, grasp the needle 2/3 of the body from the point. If the needle is held too close to the needle point, penetration qualities should be affected, if it is held too close to the hole the needle can be twisted and broken at this point. Changing the bending of the needle leads to its breakage.

Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death.

Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness.

For single patient use only - Contents supplied sterile

Do not use if the product or the pouch appears to be damaged.

USE

Inspect the packaging to confirm its integrity before opening. Open the peelable pouch and take the device.

As any suture material, the technique (knots safety), methods and associated products must be adapted by the surgeon according to the clinical context and his experience to ensure patient's safety.

Users should be familiar with surgical procedures and techniques involving synthetic non-absorbable sutures, as risk of wound dehiscence may vary with the site of application and the suture material used.

Acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

This product can be used only by or on the order of a physician.

Any device which has been used or partially used must be destroyed in proper incinerator in accordance with the local or national regulations concerning the disposal of hospital waste.

ADVERSE EVENTS

Side effects associated with the use of this medical device include:

- wound dehiscence,
- calculi formation in urinary or biliary tractus when prolonged contact with salt solutions occurs,
- Sinus formation
- Risk, if existing microbiological infection, that suture, as any foreign substance, accentuate infection,
- Minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body,
- Transitory local irritation.

Some complications may result from the surgical technical itself

VALIDITY / STORAGE CONDITIONS

Do not use after the expiry date.

Recommended storage conditions: must be stored in the original packaging, below 40°C, in a dry place, away from light.

SYMBOLS USED ON LABELLING

	= PLEDGET
	= Medical device
	= Unique Device Identifier
	= Do not reuse
	= Do not resterilize
	= Use by date (YYYY-MM-DD)
	= Date of manufacture (YYYY-MM-DD)
	= Manufacturer
	= Batch code
	= Do not use if package is damaged
	= Keep away from sunlight
	= Keep dry
	= Upper temperature limit
	= Sterilized using ethylene oxide
	= Single sterile barrier system
	= Single sterile barrier system with protective packaging inside
	= Catalogue number
	= Caution
	= Consult instructions for use
	= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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