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Circular Stapler EN

GBHEM

READ THE FOLLOWING INFORMATION THOROUGHLY BEFORE USING PRODUCT.

IMPORTANT NOTES

This booklet is NOT a reference to surgical techniques but a supporting material for the use of this product, which was designed, tested and produced merely for single patient use. It is NOT suggested that this device be reused or reprocessed or resterilized due to the risks of failure, patient injury, contamination as well as patient infection.

DESCRIPTION

The GOLD BRIDGE GBHEM Hemorrhoidal Circular Stapler places circular, double staggered rows of titanium staples in the tissue and resects the excess tissue, thus creating a circular anastomosis. The GBHEM set includes a hemorrhoidal circular stapler, Suture Threader, Circular Anal Dilator and Purse-string Suture Anoscope. There are 2 stapler diameters in common use, 32 mm, or 34 mm.

INDICATIONS

The GOLD BRIDGE GBHEM Hemorrhoidal Circular Stapler and accessories have application throughout the anal canal to perform surgical treatment of hemorrhoidal disease.

CONTRAINDICATIONS

1. Do not use the instrument if the combined compressed tissue thickness is greater than 1.5 mm, or the internal diameter of the rectum does not accommodate the instrument and accessories. If the instrument is used on tissue greater than 1.5 mm in thickness, an inadequate mucosal repair and inadequate hemostasis could result.
2. Do not use the instrument on ischemic or necrotic tissue.

WARNINGS AND PRECAUTIONS

1. Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
2. Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of procedure.
3. Do not attempt to release the safety until the instrument is ready to be fired.
4. Do not immerse the GOLD BRIDGE GBHEM Hemorrhoidal Circular Stapler and accessories in alcohol or any quaternary ammonium solutions.
5. Always inspect the staple line for hemostasis. Metal clips, staples, or sutures contained in the area to be stapled may affect the integrity of the stapled mucosal repair. Corrective action, if required, may include the use of sutures or electrocautery.
6. Ensure that the tissue thickness is within the indicated range, and that it is evenly distributed in the instrument. Excess tissue on one side may result in unacceptable staple

formation and/or incomplete knife cut. This may result in intraoperative leaks.

7. BEFORE FIRING, ENSURE THAT THE RED LINE IS IN THE MIDDLE OF THE GREEN RANGE IN THE INDICATOR WINDOW.

8. Ensure that the firing handle is fully squeezed to guarantee proper staple formation and cutting of tissue. **Partial or incomplete handle squeezes may result in unacceptable staple formation and/or incomplete knife cut.** This may result in intraoperative leaks.


9. Keep the staple line at least 2 cm above the dentate line.

10. Avoid excessive dilation.

11. Avoid incorporation of the underlying musculature in the resection or staple line.

12. Instruments or devices which come into contact with bodily fluids may require special disposal to prevent biological contamination.

13. Dispose of all opened instruments whether used or unused. Do Not resterilise the GBHEM Hemorrhoidal Circular Stapler or accessories. Resterilisation may compromise the integrity of the instrument or accessories, which may result in wound leakage or disruption.

A-GBHEM Hemorrhoidal Circular Stapler

B -- Circular Anal Dilator (CAD)

C -- Purse-String Suture-Anoscope (PSA)

D -- Circular Anal Dilator-Obturator

E -- Suture Threader

INSTRUCTIONS FOR USE

1. After an accurate lubrication of the anal margin the Circular Anal Dilator (CAD), with the obturator in place, is introduced with the aid of short rotary movements (clockwise and anti-clockwise). The introduction of the CAD causes the reduction of the prolapse, the anoderm and parts of the rectal mucosa. Sometimes introducing a surgical sponge into the anus before the CAD can improve the reduction of the excess anoderm. After removing the obturator, the prolapsed mucosa falls into the lumen of the CAD. As it is transparent, the CAD allows visualization of the dentate line, allowing the surgeon to check the correct positioning. The CAD is then stitched to the perineum with 4 stitches, at the cardinal points.

Key Points

- 1) The four anchoring stitches to secure the CAD can be placed out before inserting the CAD. By applying traction on the threads the insertion can be further facilitated.
- 2) Insert the well-lubricated CAD without stretching the sphincter. Apply counter-traction to

facilitate insertion.

3) It is possible to insert the obturator first, then extract it and insert the CAD with the obturator in place. This facilitates complete dilatation of the anus.

4) Reduce the external component as much as possible manually, even after the removal of the obturator, when the CAD is already in place.

5) Abort the procedure if stenosis precludes passage.

2. The Purse-string Suture Anoscope (PSA) is introduced through the CAD. This instrument will move the mucous prolapse along the rectal walls along a 270° circumference, while the mucosa that protrudes through the PSA window can be easily contained in a suture that includes only mucosa and submucosa. This suture must be carried out at least 2-3 cm above the apex of the hemorrhoids, the distance to be increased in proportion to the extension of the rectal prolapse. By rotating the PSA, it will be possible to complete a purse-string suture around the entire anal circumference.

Key Points

1) Refer to the apex of the hemorrhoids.

2) Carry out the Purse-string suture via a 2-0 monofilament on 25-30 mm curved needle.

3) At each stitch, extract the PSA then rotate it and insert it again. Rotating the PSA while fully inserted in the CAD can twist the mucosa and result in an improper asymmetric purse-string suture.

4) Do not tightly close the purse-string suture at this time. Confirm adequacy of purse-string, ensuring no gaps.

3. The GBHEM Hemorrhoidal Circular Stapler is opened to its maximum position. Its anvil is introduced beyond the purse-string. The GBHEM is slightly withdrawn; enough to ensure the purse-string can be visualized. The purse-string is then tied with a single closing knot. With the help of the Suture Threader (ST), the ends of the suture are pulled through and out of the lateral holes of the GBHEM.

Key Point

1) Always open the circular stapler to the maximum position.

2) Should the insertion of the anvil beyond the purse-string be difficult, do not use force: loosen the purse-string first, and then insert the anvil again.

3) Secure the purse-string under direct visualization.

4) To pull-out each suture thread, just place it opposite to the exit side (i.e. Suture Threader in the left hole, suture thread extended to the right side).

4. The ends of the suture are knotted externally or held by a clamp. The casing of the GBHEM is introduced into the anal canal. During the introduction, it is advisable to partially tighten the stapler.

Key Point

During this step the GBHEM should be gently pushed in, whilst the threads are pulled by the surgeon so that the prolapsed mucosa begins to be accommodated in the casing.

5. Via a moderate traction on the purse-string, the prolapsed mucosa is accommodated by the casing of the GBHEM. The instrument is then tightened to the end, by fully rotating the knob of the stapler clockwise. At this time the red line must be in the middle of the green

range. The stapler is then fired. Keeping the GBHEM in the closed position for 30 seconds after firing may act as a tamponade and promote hemostasis. The GBHEM is then slightly opened (one-half to three quarters of a revolution) and extracted. Finally, the staple line is examined using the PSA.

Key Point

- 1) At the end of the closure, the 4 cm mark should be at the level of the anal verge.
- 2) Align the stapler along the axis of the anal canal, and close the stapler while maintaining moderate tension on the purse-string.
- 3) If the patient is a woman, check the posterior vaginal wall to be certain that it has not been incorporated in the staple line.
- 4) One turn of the knob should be sufficient to open the stapler. Additional turns could cause the interposition of mucosa between the anvil and the upper edge of the GBHEM, thus requiring the extraction of the CAD and GBHEM simultaneously.
6. Inspect the staple line for bleeding and if necessary reinforce it with haemostatic stitches using 3-0 absorbable suture on a small needle. In order to ensure the ease of checks for post-op bleeding, providing adequate comfort for the patient, a long wet swab (25cm), with a suture anchored at the lower end, is completely inserted above the staple line. The swab is normally removed after 4-5 hours. Inspect the specimen to confirm that the technique has been properly performed. Accessory procedures (e.g. tag or papilla excision) can be performed either before or after stapling.

Key Points

- 1) After firing, leaving the CAD in place significantly facilitates the inspection of the suture line.
- 2)

Electrocoagulation must be avoided for the hemostasis because of the staples.

STAPLE SPECIFICATION CHART

Q) SIZE OF STAPLE

R) APPROXIMATE CLOSED STAPLE HEIGHT (HEIGHT OF "B")

S) SHELL OUTSIDE DIAMETER

T) NUMBER OF STAPLES

U) TISSUE COMPRESSION REQUIREMENT

	Q	R	S	T	U
32	 3.8mm	 1.5mm	32mm	32	1.5mm
34	 3.8mm	 1.5mm	34mm	34	1.5mm

STORE AT ROOM TEMPERATURE.

AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

DO NOT EXPOSE TO TEMPERATURES ABOVE 130° F (54° C).

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GOLD BRIDGE MEDICAL TECHNOLOGIES LTD
Shanghai International Holding Corporation GmbH