Linear Stapler  
Reloading Unit For Linear Stapler  
Agrafeuse Linéaire  
Unité de Recharge pour l'Agrafeuse Linéaire  
Lineares Klammernahtinstrument  
Ladeinheit für lineares Klammernahtinstrument

**GBLS/GBLR**

**Linear Stapler**  
**Reloading Unit For Linear Stapler**  
**EN**

**GBLS/GBLR**

**READ THE FOLLOWING INFORMATION THOROUGHLY BEFORE USING THIS PRODUCT.**

**IMPORTANT NOTES**

This booklet is NOT a reference to surgical techniques but a supporting material for the use of this product, which was designed, tested and produced merely for single patient use. It is NOT suggested that this device be reused or reprocessed or resterilized due to the risks of failure, patient injury, contamination as well as patient infection.

**DESCRIPTION**


The GOLD BRIDGE GBLS Surgical stapler places a double or triple staggered row of titanium staples and is available in 30 mm, 45 mm, 60 mm and 90 mm staple line lengths for use in different applications. Three staple sizes (2.5mm, 3.5 mm and 4.8 mm) are available to accommodate different tissues thicknesses. Each GBLS Surgical stapler can be reloaded 7 times for a total 8 firing on one patient.

**INDICATIONS**

The GOLD BRIDGE GBLS Surgical stapler can be applied in abdominal, thoracic and pediatric surgical procedures for transection or resection of tissue.

**CONTRAINDICATIONS**

1. This instrument is not designed, sold, or intended for use except as indicated.
2. The GOLD BRIDGE GBLS surgical stapler should NOT be used on tissues that would not tolerate conventional suture materials or conventional closure techniques if so judged by the surgeon,.
3. The GOLD BRIDGE GBLS surgical stapler should NOT be used on tissue that can not properly compress, or that compresses to less than the specified compression requirements. Refer to "Tissue Compression Requirement" section. If the instructions are not followed, closure failure, tissue trauma, dehiscence, tissue tearing, and displacement may occur, and/or hemostasis may not be obtained.
4. The GBLS Surgical staplers should not be used to staple liver, spleen, or similar tissue, as it may be destructive to compress these tissues.

## STAPLE SPECIFICATION CHART

- A) INSTRUMENT CODES
- B) RELOADING UNIT CODES
- C) STAPLE LINE LENGTH
- D) COLORS

A	B	C	D
GBLS90	GBLR90B	90mm	BLUE
	GBLR90G		GREEN
GBLS60	GBLR60B	60mm	BLUE
	GBLR60G		GREEN
GBLS45	GBLR45B	45mm	BLUE
	GBLR45G		GREEN
GBLS30	GBLR30B	30mm	BLUE
	GBLR30G		GREEN

- E) RELOADING UNIT CODES
- F) CLOSED STAPLE SIZE
- G) CLOSED STAPLE SIZE
- H) TISSUE COMPRESSION REQUIREMENT


E	F	G	H
GBLR90G	 4.8mm	 2.0mm	2.0mm
GBLR60G			
GBLR45G			
GBLR30G			
GBLR90B	 3.5mm	 1.5mm	1.5mm
GBLR60B			
GBLR45B			
GBLR30B			

**NOTE:** “Tissue Compression Requirement” refers to the compression requirement for each staple size. If tissue cannot be comfortably compressed to this requirement, or compressed to less than this requirement, the tissue may be too thick or thin for the staple size.

### INSTRUCTIONS FOR USE

THE GOLD BRIDGE GBLs Surgical staplers are all supplied sterile, and when supplied, the jaws are in the open position with one reloading unit in place.

1. Slip the jaws around the tissue to be transected or resected.
2. Before closing or firing the instrument, make certain that the retaining pin has been pushed forward and is properly seated, meanwhile the indicator window should present a white color. The pin will move when the white knob is manually pushed forward. If the retaining pin is not placed in position manually before closing the instrument, it will automatically seat when the instrument jaws are closed.

**CAUTION:** Proper seating of the retaining pin in the anvil hole should be confirmed visually and by feeling each side of the anvil to insure that the retaining pin is not misplaced. Firing an instrument with a misplaced retaining pin may result in improperly formed staples, which could result in leakage or disruption of the staple line.

**NOTE:** If the retaining pin is manually seated and the release button is activated prior to a complete handle squeeze, the pin will need to be manually retracted. In all other situations the pin will retract automatically.

3. Squeeze the handle to begin clamping the tissue. There is a tactile “pre-clamp” position built into the instrument where the handle can be released without returning to its original position to allow for final positioning of the tissue within the jaws, and the indicator window presents a yellow color. After the tissue is properly positioned, continue the handle stroke until the instrument is clamped. On release, the handle will return to the original position, and the indicator window presents a green color.

**CAUTION:** When the stapler is fully closed, the hand should leave the firing handle so that the firing handle can completely return to its original position. If the handle does not go back fully it may not be possible to fire the instrument or the handle may not lock after firing.

4. To fire the instrument, squeeze the handle again until it reaches a solid stop and locks in the back position when the indicator window becomes red. When fully fired, the handle will remain in the locked position until the black release button is pushed and the instrument is opened.

**CAUTION:** If the instrument handle is partially squeezed during firing then released, bleeding may occur as the instrument was not fully fired. The black release button may be used to open the instrument at any point during approximation.

5. Prior to opening the fully fired instrument, the cartridge edge may be used as a cutting guide during tissue transection.

6. Press the black release button in the rear of the instrument to open the jaws and unlock the handle. The alignment pin will automatically retract irrespective of whether it was pushed forward manually or automatically.

**NOTE:** After removing the instrument, always inspect the staple line for hemostasis. Minor bleeding can be controlled by manual sutures or electrocautery.


### **RELOADING**

1. Ensure that the instrument jaws are in the fully open position, and the indicator window presents a white color, and the alignment pin is fully retracted (the retaining pin knob must be in the back position) before the reloading unit is removed.

2. To remove the Reloading Unit, grasp the finger pads at the top of the Reloading Unit and pull it straight upward from the jaws.

Prior to reloading the instrument, any tissue and/or formed but unused staples must be cleared out from the anvil surface to safeguard proper staple formation in subsequent firings.

3. To reload the instrument, grasp the unused new Reloading Unit over the finger pads with the staple holes facing the instrument anvil. Insert the Reloading Unit into the metal cartridge housing and push firmly downward until the reloading unit clicks into position.

**NOTE:** The instrument jaws will not close in these positions: 1. the Reloading Unit has not been loaded into the jaws; 2. a fired or partially fired Reloading Unit is still in the jaw. Always make certain that a Reloading Unit of the proper size is used according to the tissue thickness. An empty Reloading Unit can be identified by the presence of white pushers in the front of the staple slots.

4. Repeat the loading procedure as needed.

**NOTE:** Each GBS Surgical stapler can be reloaded up to 7 times for a total of 8 firings on one patient.

**STORE AT ROOM TEMPERATURE.**

**AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.**

**DO NOT EXPOSE TO TEMPERATURES ABOVE 130° F (54° C).**

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**GOLD BRIDGE MEDICAL TECHNOLOGIES LTD**
Shanghai International Holding Corporation GmbH