**Endo Linear Cutter Stapler
Reloading unit for Endo Linear Cutter
Stapler**

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GBELCAS/GBELCR/GBELCCR

Endo Linear Cutter Stapler
Reloading unit for Endo Linear Cutter Stapler

EN

GBELCAS/GBELCR/GBELCCR

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION

**THOROUGHLY
IMPORTANT**

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Do not reuse, reprocess or resterilize this device.

DESCRIPTION

The Endo Linear Cutter Stapler place 2 x 3 rows of staggered titanium staples and cut between the 2 x 3 staggered rows of staples. The size of the staples is determined by the selection of the 2.5mm, 3.5mm, 4.2mm, 4.8mm or 5.0mm Reloading Unit for Endo Linear Cutter Stapler (RU). The Endo Linear Cutter Stapler will accommodate any of the reloading unit sizes that are available in the 30mm, 45mm, and 60mm lines.

The Endo Linear Cutter Stapler with either the 2.5mm, 3.5mm ,4.2mm, reloading unit is designed for introduction and use through a 12mm trocar sleeve, or larger, with the use of a converter.

When using the Endo Linear Cutter Stapler with the 4.8mm, 5.0mm reloading unit, it must be inserted into a 15mm trocar sleeve.

The instrument may be reloaded and fired up to 25 times in a single procedure.

GBELCCR series are articulating reloading units with Curved tip.

The curved tip of GBELCCR provides better visualization to ensure that target tissues are captured in difficult-to-reach area, which greatly improves the maneuverability and visualization.

NOTE: Each instrument can accommodate the 30-2.5, 3.5, 4.2, 4.8, 5.0/ 45-2.5, 3.5, 4.2, 4.8, 5.0/60-2.5, 3.5, 4.2, 4.8, 5.0 mm RU's.

INDICATIONS

The Endo Linear Cutter Stapler has applications in abdominal, gynecologic, pediatric and thoracic surgery for resection, transaction and creation of anastomosis.

CONTRAINDICATIONS

1. Do not use the 2.5mm staples on any tissue that compresses to less than 1.0mm in thickness, on any tissue that cannot comfortably compress to 1.0mm or on aorta.
2. Do not use the 3.5mm staples on any tissue that compresses to less than 1.5mm in thickness, on any tissue that cannot comfortably compress to 1.5mm or on aorta.
3. Do not use the 4.2mm staples on any tissue that compresses to less than 1.8mm in thickness, on any tissue that cannot comfortably compress to 1.8mm or on aorta.
4. Do not use the 4.8mm staples on any tissue that compresses to less than 2.0mm in thickness, on any tissue that cannot comfortably compress to 2.0mm or on aorta.
5. Do not use the 5.0mm staples on any tissue that compresses to less than 2.35mm in thickness, on any tissue that cannot comfortably compress to 2.35mm or on aorta.
6. The Endo Linear Cutter Staple instrument should not be used on tissue such as liver or spleen where compressibility is such that closure of the instrument would be destructive.
7. Do not use the Endo Linear Cutter Staple where adequacy of hemostasis cannot be verified visually after applications.
8. These devices are provided STERILE and are intended for use a SINGLE procedure only.

DISCARD AFTER USE. DO NOT RESTERILIZE.

WARNINGS AND PRECAUTIONS


1. Preoperative radiotherapy may result in changes to tissue.

These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple size. Careful consideration should be given to any pre-surgical treatment the patient may have undergone and in corresponding selection of staple size.

2. Always include the combined thickness of the tissue and of any staple line reinforcement material use when choosing the proper staple cartridge.
3. When using the Endo Linear Cutter Staple with a 4.8, 5.0 reloading unit, the instrument MUST be inserted into a 15mm trocar. A smaller size trocar will not accept the 4.8, 5.0 reloading unit.
4. Always close the jaws of the Endo Linear Cutter Staple prior to introducing and removing the stapler from the trocar sleeve.
5. After firing, always inspect the staple line for hemostasis. Minor bleeding may be controlled by electrocautery or manual sutures.
6. Placement of tissue proximal to the tissue stops (on the RU) may result in stapler malfunction. Any tissue extending beyond the cut mark will not be transected.
7. When using the stapler more than once during a SINGLE surgical procedure, be sure to remove the empty RU and reload a new one. A safety interlock is provided that prevents an empty RU from being fired a second time. Do not attempt to override the safety interlock.
8. When positioning the stapler on the application site, ensure that no obstructions, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
9. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Prior to performance of any endoscopic procedures, consult the medical literature relative to techniques, complications and hazards.
10. A thorough understanding of the principles involved in laser and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and operator(s), and damage to the instrument.
11. When endoscopic instruments and accessories from different manufactures are employed together in a procedure, verify compatibility and ensure that electrical isolation or grounding is not compromised.
12. The anvil must be completely visible, (past the trocar sleeve) prior to opening the RU within the body cavity.
13. When using a staple line buttressing material, follow the instructions provided by the manufacturer of the buttress material, as performance of the stapler may be affected when using buttress materials.
14. The instrument and reloading unit are provided STERILE and are intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZE.
15. Do not attempt to load RU while squeezing ring handle.

① SCHEMATIC VIEW


- A) PIN
- B) UNLOAD/UNLOCK BUTTON
- C) SHAFT
- D) ARTICULATING LEVER
- E) ROTATION COLLAR
- F) GREEN BUTTON
- G) RETURN KNOB
- H) HANDLE
- I) RELOADING UNIT(RU)(30,45,60)
- J) LOAD ALIGNMENT INDICATORS
- K) SHIPPING WEDGE
- L) ANVIL

- M) TISSUE STOP
- N) END OF STAPLE LINE
- O) END OF CUT LINE
- P) INCREMENT MARKINGS
- Q) LOWER CLAMP BUTTON
- R) STAPLE CARTRIDGE

② LOADING

2

1. The RU is packaged in the open position. Do not attempt to close the RU.

WARNING: SELECT A RU WITH THE APPROPRIATES STAPLES SIZE FOR THE TISSUE THICKNESS.-OVERLY THICK OR THIN TISSUE MAY RESULT IN UNACCEPTABLE STAPLE FORMATION. ALWAYS INCLUDE THE COMBINED THICKNESS OF THE TISSUE AND OF ANY STAPLE LINE REINFORCEMENT MATERIAL IN USE WHEN CHOOSING THE PROPER STAPLE CARTRIDGE.

CAUTION: Do not attempt to remove the shipping wedge until the RU is loaded into the instrument.

2. Ensure that the knobs on the instrument are pulled back completely and the articulation arm is neutral to the instrument.

G) RETURN KNOBS

3. To load the Endo Linear Cutter Stapler with appropriate RU, insert the pin located at the distal end of the instrument shaft into the RU. Ensure that the LOAD alignment indicator on the RU aligns with the LOAD alignment indicator on shaft. Push the RU in and twist clockwise 45° relative to the instrument, so that the RU will lock into place. The LOAD alignment indicator on the instrument shaft will align with the LOAD alignment indicator on the RU.

J1) LOAD ALIGNMENT INDICATOR (RU)

J2) LOAD ALIGNMENT INDICATOR (SHAFT)


4. Remove the shipping wedge from the RU prior to inserting the instrument into the trocar.

CAUTION: Do not clamp instrument prior to removing shipping wedge.

5. To confirm proper loading, cycle the instrument after loading the RU. Squeeze the handle once to close the jaws of the RU. Pull back on the return knob and confirm that the RU jaws open fully.

③ UNLOADING

3

1. To unload a RU from the stapler, the articulating lever must be in the neutral position. Ensure that the jaws of the RU are open by pulling the return knobs back completely. Pull the UNLOAD/UNLOCK button (located on the underside of the shaft) back towards the instrument, twist the RU counterclockwise 45° and remove the RU from the shaft of the instrument.

D) ARTICULATING LEVER

B) UNLOAD/UNLOCK BUTTON


④ INSTRUCTIONS FOR USE

NOTE: The jaws of the RU must be closed prior to introducing the instrument into the trocar sleeve. To do so, squeeze the handle.

1. Insert the Endo Linear Cutter Stapler into an appropriately sized trocar sleeve, or larger, with the use of a converter.

CAUTION: The anvil must be completely visible, (past the trocar sleeve) prior to opening the RU within the body cavity.

The instrument shaft can be rotated 360° by using the rotation collar. The articulated RSU's of GBELCAS series can be articulated at 22° and 45°, in both directions, by using the articulating lever.

NOTE: When using the Endo Linear Cutter Stapler instrument with a 4.8 or 5.0 RU, the instrument **MUST** be inserted into a 15mm trocar. A smaller size trocar will not accept the 4.8 RU.

2. Once inside the body cavity, open the jaws of the instrument by pulling the return knobs completely back.

CAUTION: Do not squeeze the instrument handle while pulling back the knobs.

3. Apply the Endo Linear Cutter Stapler across the tissue to be transected. Caution: Ensure that no obstructions (such as clips) are incorporated in the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

The instrument will not cut tissue beyond the black cut mark indicated on the RU. More than one application of the Endo Linear Cutter Stapler may be necessary for tissue exceeding the length of the RU (30mm, 45mm, or 60mm).

S) CUTMARK

CAUTION: Placement of tissue proximal to the tissue stops (on the RU) may result in stapler malfunction. Any tissue extending beyond the cut mark will not be transected.

4. Close the jaws of the instrument across the tissue to be transected by squeezing the handle completely. The stapler is equipped with a safety interlock; the instrument will not fire the staples and cut tissue unless the green button is pushed.

CAUTION: A safety interlock is provided that prevents an empty RU from being fired a second time. Do not attempt to override the safety interlock.

The jaws of the instrument may be repositioned on the tissue prior to firing by pulling the return knob completely back, allowing the jaws to open.

5. In order to fire the instrument, push the green button. Squeeze the handle sequentially until the oval clamp cover reaches the distal end of the cartridge slot, and the handle locks.

Sequential squeezes of the handle are required to fully fire the RU.

The total number of squeezes is relative to the length of the RU

(30, 45 or 60).

FAILURE TO COMPLETELY FIRE THE RU WILL RESULT IN AN INCOMPLETE CUT AND/OR INCOMPLETE STAPLE FORMATION, WHICH MAY RESULT IN POOR HEMOSTASIS.

F) GREEN BUTTON

Q) LOWER CLAMP BUTTON

6. Once the instrument has been completely fired, open the jaws by pulling the return knob completely back. Gently remove the instrument from the tissue. The site should be hemostasis following removal of the instrument. Minor bleeding may be controlled using electrocautery or manual sutures.

7. After completely firing the RU, close jaw of the instrument and remove the Endo Linear Cutter Stapler from the body cavity to unload the RU from the instrument.

NOTE: Do not attempt to insert or remove the instrument from the trocar sleeve if the instrument is in the articulated position.

The Endo Linear Cutter Stapler may be reloaded and fired up to 25 times in a single procedure.

⑤ STAPLE SPECIFICATIONS

T) INSTRUMENT REORDER CODES

U) RU REORDER CODES

V) STAPLE LINE LENGTH

W) COLOR

X) OPEN STAPLE SIZE

Y) CLOSED STAPLE SIZE

T	U	V	W	X	Y
GBELCAS 06 GBELCAS 16 GBELCAS 26	GBELCR3 0W	30mm	White	 2.5mm	1.0mm
	GBELCR3 0B	30mm	Blue	 3.5mm	1.5mm
	GBELCR3 0D	30mm	Gold	 4.2mm	1.8m06 m
	GBELCR3 0G	30mm	Green	 4.8mm	2.0mm
	GBELCR3 0N	30mm	Black	 5.0mm	2.35mm
	GBELCR4 5W	45mm	White	 2.5mm	1.0mm
	GBELCR4 5B	45mm	Blue	 3.5mm	1.5mm
	GBELCR4 5D	45mm	Gold	 4.2mm	1.8mm
	GBELCR4 5G	45mm	Green	 4.8mm	2.0mm
	GBELCR4 5N	45mm	Black	 5.0mm	2.35mm
	GBELCR6 0W	60mm	White	 2.5mm	1.0mm
	GBELCR6 0B	60mm	Blue	 3.5mm	1.5mm
	GBELCR6 0D	60mm	Gold	 4.2mm	1.8mm
	GBELCR6 0G	60mm	Green	 4.8mm	2.0mm
	GBELCR6 0N	60mm	Black	 5.0mm	2.35mm

T	U	V	W	X	Y
GBELCAS 06 GBELCAS 16 GBELCAS 26	GBELCCR 30W	30mm	White	 2.5mm	1.0mm
	GBELCCR 30B	30mm	Blue	 3.5mm	1.5mm
	GBELCCR 30D	30mm	Gold	 4.2mm	1.8m06 m
	GBELCCR 30G	30mm	Green	 4.8mm	2.0mm
	GBELCCR 30N	30mm	Black	 5.0mm	2.35mm
	GBELCCR 45W	45mm	White	 2.5mm	1.0mm
	GBELCCR 45B	45mm	Blue	 3.5mm	1.5mm
	GBELCCR 45D	45mm	Gold	 4.2mm	1.8mm
	GBELCCR 45G	45mm	Green	 4.8mm	2.0mm
	GBELCCR 45N	45mm	Black	 5.0mm	2.35mm
	GBELCCR 60W	60mm	White	 2.5mm	1.0mm
	GBELCCR 60B	60mm	Blue	 3.5mm	1.5mm
	GBELCCR 60D	60mm	Gold	 4.2mm	1.8mm
	GBELCCR 60G	60mm	Green	 4.8mm	2.0mm
	GBELCCR 60N	60mm	Black	 5.0mm	2.35mm

STORE AT ROOM TEMPERATURE.

**AVOID PROLONGED EXPOSURE TO ELEVATED
TEMPERATURES.**

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GOLD BRIDGE MEDICAL TECHNOLOGIES LTD
Shanghai International Holding Corporation GmbH