INSTRUCTIONS FOR USE POLYMESH PLUG

PRODUCT DESCRIPTION:

Polymesh Plug is a pre-formed, three-dimensional device constructed of a fluted outer layer of mesh and multiple inner layers of mesh attached at the tip and a separate flat, pre-shaped onlay mesh which packaged together.

INDICATIONS:

Polymesh Plug is indicated for the repair of groin hernia defects.

CONTRAINDICATIONS:

Do not use Polymesh Plug in infants or children, whereby future growth will be compromised by use of such mesh material.

WARNINGS:

- 1.The use of any permanent mesh or patch in a contaminated or infected wound could lead to fistula formation and/or extrusion of the mesh.
- 2.If an infection develops, treat the infection aggressively/Consideration should be given regarding the need to remove the mesh. An unresolved infection may require removal of the device.
- 3.Polymesh Plug is supplied sterile. Inspect the packaging to be sure it is intact and undamaged prior to use.
- 4.Polymesh Plug is for single use only. Do not resterilize or reuse any portion of the Polymesh Plug
- 5.To avoid injury, careful attention is required if fixating the mesh in the presence of nerves or vessels.
- 6.Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient. Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection including but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient or end user.

PRECAUTIONS:

- 1.Only physicians qualified in the appropriate surgical techniques should use this prosthesis.
- 2.Please read all instructions prior to use.
- 3.Care should be taken to ensure that the mesh is adequately fixated to the uncompromised tissue of the inguinal floor. If necessary, additional fasteners and/or sutures should be used.
- 4.Monofilament sutures are recommended to properly secure Polymesh Plug.

DIRECTIONS FOR USE:

For indirect hernias, Polymesh Plug is inserted through the internal ring, tapered end first, and positioned just beneath the crura.

Polymesh Plug should be sutured in place with interrupted sutures. If the overall bulk of Polymesh Plug is deemed excessive, then some of the inside "petals" may be trimmed. For direct hernias, the defect is circumscribed at its base and the contents fully reduced. Polymesh Plug is then inserted through the defect and routinely secured in position with interrupted sutures. The placement of two or more plugs has been used successfully for large direct and pantaloon hernias. In such a situation where plugs abut, they should be joined together with sutures. In femoral hernias a small Polymesh Plug, with the interior petals removed, is secured in position with interrupted sutures after the sac has been reduced. Once the sac is reduced, the plug is placed such that the top of the plug lies flush with the opening of the femoral canal. Polymesh Plug is secured to the surrounding tissue with sutures. Care must be taken when suturing the plug to prevent trauma to the femoral vessels and nerve. An onlay patch is used only if there is adequate space for it to be placed. In all groin hernias, the flat, pre-shaped onlay is used for reinforcement. It is placed without sutures on the anterior surface of the posterior wall of the inguinal canal. The tails of the onlay patch will help reinforce the internal ring. If desired, the tails can be brought together with sutures. Nonabsorbable, monofilament sutures are recommended. Fibroblastic proliferation and collagen formation will penetrate the mesh and secure it in place.


ADVERSE REACTIONS:

Possible complications include seroma, adhesions, hematoma, inflammation, extrusion, fistula formation and recurrence of the hernia or soft tissue defect.

TRACEABILITY:

A traceability label is attached to every prosthesis package which identifies the type, size, and lot number of the prosthesis. This label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted.

APPLICATION:

Symbols Used on Labelling

	Manufacturer		Caution
	Date Of Manufacture		Temperature Limit
	Use-By Date		Humidity Limit
	Size/Model		Do Not Re-Use
	Product Code		Do Not Resterilize
	Batch Number		Do Not Use If Package Is Damaged
	Sterilized Using Ethylene Oxide		Keep away from direct sunlight
	Consult Instructions For Use		Keep Dry
	Notified Body Number		Authorized representative in the European Community

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Betatech Medical Cihazlar San. Müm. İç ve Dış Tic. Ltd. Şti.