INSTRUCTIONS FOR USE

REDON DRAINAGE TUBE

Materials Used

- Poly Vinyl Chloride (PVC).

Indications

- Collection of fluids from surgical wound under negative pressure.

Contraindications

- Fistulas to organs or body cavities.
- Necrotic tissue in escher bacteria.
- Osteomyelitis (untreated).
- Malignancy in the wound.

Directions for Use

- Place the perforated tube in the wound using suitable method e.g. with the help of curved needle. Ensure that the entire perforated part of the catheter is well within the wound to avoid any air leakage during drainage.
- The length of the perforation should be reduced according to the size of the wound.
- Catheter is marked 5 cm from the last perforation to ascertain the depth of catheter in the wound.
- Connect the tube with a suitable suctioning tube for drainage.

Cautions

- Do not use if package is opened or damaged.
- The product should be used only by qualified healthcare professionals.
- Do not re-sterilize. Discard the set after single use.
- Reuse and cleaning of product may alter their structural and mechanical properties. It may lead to infection or other illness/injury.
- Store in between 5°C to 35°C, avoid excessive heat, protect from direct sunlight and moisture.
- Any device used must be disposed of in an appropriate incinerator in accordance with local or national regulations concerning the disposal of hospital waste.
- Do not tighten excessively when fixing the drain to the skin. Caution should be taken not to perforate the drain with the needle or with other sharp instrument what would result in damaging the tube. Do not leave in the patient longer than necessary.
- Poly Medicure Limited will not be responsible for any direct incidental or consequential damages resulting from reuse of product.

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POLY MEDICURE LIMITED
 Plot No. 104-105, Sector 59, HSIIDC Indl. Area,
 Ballabgarh, Faridabad-121004 (Haryana), INDIA

EC REP Obelis s.a.,
 Boulevard Général Wahis 53,
 1030, Brussels, BELGIUM.
 Tel : +(32) 2. 732.59.54
 E-mail : mail@obelis.net

- Cautions
- Consult Instructions for use
- Catalogue reference
- For single use only
- Do not Resterilize
- Do not use if packaging or product has been damaged or contaminated
- Batch Number
- Date of Manufacturing
- Use by / Expiry date
- Sterilised by Ethylene Oxide
- Single Sterile Barrier System with Protective Packaging Inside
- European Authorised Representative
- Manufactured by
- Pyrogen Free
- DEHP FREE
- LATEX FREE
- Temperature Limit
- Keep Away from Sunlight
- Keep Dry
- Medical device
- Entity distributing the Medical Devices into EU
- CE marking: The product is in conformance with the basic requirements of European Directive 93/42/EEC concerning Medical Devices

Rev. 03-21/08

Note : IFU will be put inside the duplex box

| | | | | |
|--|----------------------------------|---------------------|----------------------|---------------------|
| Doc. No. Peters/001/002-I | Date : 04-02-20 | Design by | Checked by | Approved by |
| IFU Design for Redon Drainage Tube Front Side | Size : 180x130mm | Designer | QA Person | QA Incharge |
| | P. Area : | | | |
| Color Break | Specifications : | Issue No. 04 | Scale : 1:1 | Sheet 1 of 1 |
| Black | 1. 2. | | | |
| Rev. | Amendment | Remarks | Date of Amdt. | |
| 1. | Updated new CE 0123 | | | 21-05-20 |
| 2. | MD symbol added | | | 21-06-21 |
| 3. | Added Peter's new logo & address | | | 27-08-21 |
| 4. | | | | |

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