Suction Catheter

Endo Bronchial Suction Catheter

CE 0123

	Materials	References	Radio-opacity X.R.O Line	Characteristics
Suction Catheter	PVC	06S020xx 06S021xx	No	xx (depends on references): 6, 8, 10, 12, 14, 16 or 18 Lengths: 30 cm or 50 cm Terminal orifice with lateral eye Straight No mark

DESCRIPTION

PVC suction catheter

INDICATION

To clean the secretion from the respiratory tract

CONTRAINDICATIONS

Not to be used in patients with known hypersensitivity to any of the material used

PRECAUTIONS

Do not use if package is damaged. Destroy after use. The design of this device doesn't allow its reuse or a new sterilization.

These devices should only be used by fully qualified personnel.

Technical procedure:

Connect the catheter to a vacuum system with a depression corresponding to 150 cm of water.

Open the cap of the swivel connector.

Insert the distal end of the catheter into the cannula or intubation catheter.

Proceed to the secretions aspiration

Discard the catheter and rinse the suction tube after use.

Store in between 5°C to 35°C, avoid excessive heat, protect from direct sunlight and moisture.

ADVERSE EVENTS

None adverse event known.

SHELF LIFE / PREVENTION / ELIMINATION

Recommended storage conditions: away from moisture and direct heat.

Do not use after expiration date.

Any device used or started must be disposed of an appropriate incinerator in accordance with local or national regulations concerning the disposal of hospital waste.

SYMBOLS USED ON LABELLING

	Consult instructions for use
	Date of manufacture
	Do not reuse
	Use by date
	Batch code
	Manufactured by
	Sterilized using ethylene oxide
	Single Sterile Barrier System
	Do not Resterilize
	Pyrogen Free
	Product reference
	Caution, consult instructions for use
	Do not use if package is damaged
	Temperature Limit
	Keep Away from Sunlight
	Keep Dry
	Medical Device
	Entity distributing the Medical Devices into EU

CE 0123

CE mark, products comply with the essential requirements of the 93/42/CEE European Directive relative to medical devices.

Rev. 02-21/08

POLY MEDICURE LTD.
Plot No. 104-105, Sector 59, HSIIDC Indl. Area
Faridabad-121004 (Haryana), INDIA

Peters
SURGICAL
Distribué par / Distributed by :
PETERS SURGICAL
Immeuble AURELIUM - 1 cours de l'Île Seguin
92100 BOULOGNE-BILLANCOURT
FRANCE
www.peters-surgical.com

EC REP Obelis s.a.,
Boulevard Général Wakis 53,
1030, Brussels, BELGIUM,
Tel : +(32) 2. 732.59.54
E-mail : mail@obelis.net

Doc. No. Peters Surgical	Date : 11-02-20	Design by	Checked by	Approved by
IFU Design for Suction Catheter 30cm & 50cm - Front Side	Size : 210x297mm P. Area :			
Color Break Black	Specifications : 1. 2.	Designer	QA Person	QA Incharge
Rev.	Amendment	Issue No. 03	Scale : 1:1	Sheet 1 of 1
1.	Updated new CE 0123			21-05-20
2.	Added Peter's new logo & address and MD symbol			27-08-21
3.				
4.				

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