Colle chirurgicale synthétique IFABOND®

Français (French - fr)	1
English (en)	3
Shqip (Albanian - sq)	5
العربية (Arabic - ar)	7
Čeština (Czech - cs)	9
Hrvatski (Croatian - hr)	11
Dansk (Danish - da)	13
Suomi (Finnish - fi)	15
Deutsch (German - de)	17
Ελληνικά (Greek - el)	19
Magyar (Hungarian - hu)	21
Italiano (Italian - it)	23
Latviešu (Latvian - lv)	25
Lietuvių K. (Lithuanian - lt)	27
Македонски Јазик (Macedonian - mk)	29
Português (Portuguese - pt)	33
Română (Romanian - ro)	35
Русский (Russian - ru)	37
Español (Spanish - es)	43
Svenska (Swedish - sv)	45
Türkçe (Turkish - tr)	47

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Synthetic surgical glue 0459

DEVICE IDENTIFICATION

IFABOND® is a sterile absorbable class III implantable medical device. IFABOND® is packaged in a hermetically sealed polypropylene vial. Each vial contains 0.5 ml, 1 ml or 1.5 ml of glue and is packed in a double sterile blister supplied with a syringe and a single-use needle for extracting the glue from the vial.

INTENDED USE

IFABOND® is intended to be used during surgical procedures for its adhesive action and its role as hemostatic sealant agent. IFABOND® is intended for use in open and laparoscopic surgery.

INDICATIONS

IFABOND® is indicated for its adhesive action and its role as hemostatic sealant agent for:

- Mesh fixation in digestive surgery for treatment of hernia and in uro-gynecological surgery for sacrocolpopexy.
- Staple-line reinforcement in bariatric surgery for sleeve and closure of mesenteric defects in bariatric surgery for bypass.

CONTRAINDICATIONS

- The glue should not be applied on fragile and/or damaged tissue that may not withstand thermal, inflammatory or mechanical stress.
- The use of glue is not suitable for ophthalmic surgery with PETERS SURGICAL's existing applicators.
- The glue should not be applied on cerebral tissue.
- Use of the glue is contra-indicated in orthopedic surgery.
- The glue should not be applied to incisions or openings with significant tension without being combined with other fixation or closure methods.
- Use of the glue is contra-indicated for any infected wound or where there is a risk of gangrene.
- The glue should not be used on patients with a known hypersensitivity to cyanoacrylate and/or its related degradation products.
- The glue should not be used on pregnant women, or on patients with a systemic pre-operative infection, with uncontrolled diabetes, or with conditions known to interfere with the healing process.

TARGET POPULATION

IFABOND® can only be used by a surgeon in surgical techniques and procedures involving the use of surgical glue. The procedures claimed with the device have been demonstrated in an adult population.

COMPOSITION

IFABOND® is composed of:

- ≥ 98% n-hexyl cyanoacrylate monomer,
- <2% methyl cyanoacrylate monomer and anionic stabiliser
- 0.1% free-radical stabilizer

After application, the anionic stabiliser is eliminated in its gaseous state.

MODE OF ACTION

The mode of action of IFABOND® is based on the chain polymerization of n-hexyl cyanoacrylate in contact with tissue and blood, creating strong mechanical bonds to the tissue and a homogeneous film acting as a hemostatic barrier. During application, the liquid glue polymerizes in approximately a minute in contact with the tissue. After polymerization, the glue has no longer any adhesive power with the adjacent tissue. This is followed by a progressive encapsulation of the glue by the fibrous connective tissue and a progressive resorption of the medical device. The hydrolytic and enzymatic absorption process leads to a degradation of the glue into cyanoacetic acids, alcohols and formaldehyde, all of which are metabolized by the body. The duration of resorption varies according to the type of tissue and the quantity of glue applied.

CLINICAL PERFORMANCES AND CLINICAL BENEFITS

Until the EUDAMED database is fully operational, the summary of safety and clinical performance is available on request.

The clinical performance for IFABOND® claimed by the manufacturer is adhesion between tissues or between mesh and tissue.

The clinical benefits of IFABOND® intended by the manufacturer are the following:

- Reduction of operative time.
- Reduction of hospital stay.
- Reduction of reoperation rate for digestive surgery and bariatric surgery.
- Improvement of the quality of life for digestive and uro-gynecological surgery.

INSTRUCTIONS FOR USE

Preparation of IFABOND® glue:

- Check the integrity of the packaging. Do not use if package is damaged.
- Take out the blister from refrigerator until it reaches room temperature before use.
- Open the external blister and place the sterile internal blister directly on the operating table.
- Check the fluidity and transparency of the glue.
- Open the internal blister, take out the sterile vial,
- Take out the needle and the syringe to collect the glue,
- Do not puncture the cap of the vial with the needle.

Application of IFABOND® glue :

- Unscrew the cap of the IFABOND® vial and collect the product using the syringe and needle supplied. Do not purge the syringe and needle after collection. If no applicator is used, apply the glue with the needle on the syringe over the concerned area.
- If an applicator is used, remove the puncture needle, screw the applicator onto the filled syringe. The glue can be delivered continuously or dropwise, depending on the type of applicator used.
- Adjust the pressure on the plunger of the syringe to control the distribution of the glue and apply it in a thin layer or dropwise to obtain maximum adhesion and sealing. Do not overlap or apply multiple drops of glue to the same spot. Incorrect application will result in increased polymerization time and decreased adhesion and sealing efficiency.
- It is recommended to wait approximately 1 min for complete polymerization to ensure adhesion and sealing.

After application or in case of obstruction of the application material, the residual volume of glue should be discarded with the packaging and accessories following the instructions defined in the STORAGE/HANDLING/ELIMINATION section. If necessary, all the material must be changed.

PRECAUTIONS FOR USE

Use:

- IFABOND® should not be in contact with polycarbonate devices or instruments.
- To ensure correct removal of the adhesive, only use the needle and syringe supplied by PETERS SURGICAL in the IFABOND® double blister pack. The needle must not come into contact with the patient's tissue.
- To ensure proper application of the glue, use only IFABOND® applicators supplied by PETERS SURGICAL. IFABOND® must not be diluted or mixed with other substances.
- Apply the glue in a continuous line or dropwise to obtain maximum adhesion and sealing. It is important not to overlap or apply multiple drops of glue to the same spot.
- Do not use IFABOND® glue without controlling the application (volume and flow). If too much is accidentally applied, remove it quickly with a wet compress within seconds of application or rinse thoroughly with saline.
- After polymerization, do not peel off IFABOND® to avoid tissue damage.
- Adhesion and sealing are weakened on tissues covered with a significant film of water or blood.
- Any product deposited in excess and polymerized can lead to the formation of small agglomerates of glue. These must be removed mechanically to prevent trauma to adjacent tissues and to avoid localized tissues inflammation.
- Avoid contact of the glue with textiles, surfaces, or instruments (latex or steel) that cannot be rinsed thoroughly

with lukewarm water or easily cleaned with a solvent such as acetone. Do not use acetone on tissues.

- In case of accidental contact with eyes: rinse thoroughly with saline solution.

User :

- The user will select the appropriate packaging (IB05; IB; IB+) for its surgical procedure. The user will also select the appropriate method of application of the device and, if applicable, the type of applicator to ensure the best clinical results.
- During preparation, the user may be exposed, through contact, to an exothermic reaction linked to polymerization.

Patients:

- The application of adhesive can lead to an exothermic reaction at the point of application due to polymerization.

UNDESIRABLE EFFECTS

Any serious incident arising in connection with the use of this device must be reported to the competent authority of the country in which the user and/or patient is established and to the manufacturer or its local representative.

The following side effects associated with the use of IFABOND® Surgical Glue are:

- Allergic reaction to a component of the medical device
- Moderate tissue inflammatory response characteristic of the foreign body response.
- Urinary bladder concretion.

Some complications or post-surgical pain may result from the surgical technique itself.

STORAGE/HANDLING/ELIMINATION

IFABOND® should be stored in its original packaging, away from direct sunlight, humidity and at a temperature not exceeding 5°C. Under the recommended storage conditions, IFABOND® has an expiration date of 2 years. Do not use the product after the expiration date.

Any used or opened device should be discarded in an appropriate incinerator in accordance with local or national regulations for the disposal of hospital waste. After use, the needle must be discarded in an appropriate container for sharps waste from care activities with infectious risks and similar.

SINGLE USE DEVICE

Device supplied for single use; contents supplied sterile.

Do not reuse, reprocess or resterilize.

Re-sterilization or reuse may:

- Compromise the structural integrity of the device and/or lead to device failure, which may result in polymerization.
- Generate contamination of the device and/or cause cross-infection of the patient, including, but not limited to, transmission of infectious diseases from one patient to another.

INFORMATION TO BE PROVIDED TO THE PATIENT

The surgeon is required to inform the patient of all contraindications, side effects, risks associated with the use of this medical device.

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