**Linear Cutter Stapler
Reloading Unit For Linear Cutter Stapler**

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GBLCS/GBLCR

**Linear Cutter Stapler
Reloading Unit For Linear Cutter Stapler
EN**

GBLCS/GBLCR

**READ THE FOLLOWING INFORMATION THOROUGHLY
BEFORE USING PRODUCT.**

IMPORTANT!

This booklet is NOT a reference to surgical techniques but a supporting material for the use of this product, which was designed, tested and produced merely for single patient use. It is NOT suggested that this device be reused or reprocessed or resterilized due to the risks of failure, patient injury, contamination as well as patient infection.

DESCRIPTION

The GOLD BRIDGE GBLCS Surgical Stapler Reloading Staples Unit place two double staggered rows of titanium staples and simultaneously cut and divides tissue between the two double rows. The GBLCS Surgical Stapler is available in 60 mm, 80 mm, and 100 mm lengths, and in two staple sizes, 3.8mm and 4.8mm, to accommodate various tissue thickness.

Each GBLCS Surgical Stapler may be reloaded up to 7 times for a total of 8 firings on one patient in one operation.

INDICATIONS

The GOLD BRIDGE GBLCS Surgical Staplers can be applied in abdominal, gynecological, pediatric and thoracic surgical procedures for resection, transection and creation of anastomosis.

CONTRAINDICATIONS

1. This instrument is not designed, sold, or intended for use except as indicated.
2. The GOLD BRIDGE GBLCS surgical stapler should NOT be used on tissues that would not tolerate conventional suture materials or conventional closure techniques if judged so by the surgeon,.
3. The GOLD BRIDGE GBLCS surgical stapler should **NOT** be used on tissue that can not properly compress, or that compresses to less than the specified compression requirements. Refer to "Tissue Compression Requirement" section, if the instructions are not followed, closure failure, tissue trauma, dehiscence, tissue tearing, and displacement may occur, and/or hemostasis may not be obtained.
4. The GBLCS Surgical staplers should not be used to staple liver, spleen, or similar tissue, as it may be destructive to compress these tissues.

WARNINGS AND PRECAUTIONS

1. Preoperative radiotherapy may result in tissue thickness to exceed the indicated range for the selected staple size. Pre-surgical treatment over the patient deserves high attention. The staple size should be carefully selected.
2. The GBLCS Surgical Staplers can be reloaded up to 7 times for a total 8 firing on one patient in one operation.

Additional applications may cause improperly formed staples resulting in leakage or disruption of the staple line.

3. When using the instrument more than once during the same procedure, be sure that the anvil is clear of tissue, blood and staples following each application.

4. When used in the abdomen, make sure no omental or mesenteric vessels are caught between the cartridge and anvil assembly prior to locking the instrument.

5. The GBLCS Surgical Staplers and Staples are provided **STERILE** and are intended for use in a **SINGLE** procedure only. **DISCARD AFTER USE. DO NOT RESTERILIZE.**

6. Ensure that GBLCS Surgical Staples are selected with the appropriate staple size for the tissue thickness. Overly thick or thin tissue may result in unacceptable staple formation.

7. After firing, always inspect the staple line for hemostasis. Minor bleeding can be controlled by electrocautery or manual sutures.

8. Failure to advance the handle completely may result in incomplete staple formation, which may compromise the integrity of the staple line.

9. **Do not turn the firing knob during firing. Turning the knob during firing may cause damage and possible instrument malfunction.**

CODE	STAPLE LINE LENGTH	CUTTING LINE LENGTH	STAPLE SIZE	CLOSED STAPLE SIZE	TISSUE COMPRES- SION REQUIRE- MENT	COLOR
GBLCS60B	64	60	3.8mm	 1,5 mm	1.5mm	BLUE
GBLCS60G	64	60	4.8mm	 2,0 mm	2.0mm	GREEN
GBLCS80B	84	80	3.8mm	 1,5 mm	1.5mm	BLUE
GBLCS80G	84	80	4.8mm	 2,0 mm	2.0mm	GREEN

GBLCS100B	104	100	3.8mm	 1,5 mm	1.5mm	BLUE
GBLCS100G	104	100	4.8mm	 2,0 mm	2.0mm	GREEN

INSTRUCTIONS FOR USE

FIRING INSTRUCTIONS CAUTION: Ensure that a GBLCS is selected with the appropriate staple size for the tissue thickness. Overly thick or thin tissue may result in unacceptable staple formation.

1. Open the instrument by pushing in on the release button on the cartridge fork lever handle (located at end of instrument).

2. Remove the shipping wedge from the cartridge channel by lifting and pulling it out at the finger tabs. (Remove shipping wedge after the cartridge is fully loaded. The cartridge will float up and down in final loaded position.)

3. Position the tissue to be stapled between the instrument cartridge and anvil fork; or, if the tissue is to be anastomosed, place the cartridge fork in one lumen and the anvil fork into the other lumen as appropriate. Align tissue edges equally on cartridge and anvil. This may be accomplished with the instrument halves separated or hinged at the rear. Close the instrument by moving the lever handle towards the instrument body until the audible click is heard.

4. With the instrument securely closed, rotate the firing knob to either side of the instrument.

NOTE: In its pre-fire position, the firing knob should turn to either side of the instrument.

5. Place a thumb behind the firing knob and two fingers on the instrument shoulders. Fire the instrument by sliding the firing knob forward to a complete stop.

CAUTION: Do not turn the firing knob during firing. Turning the knob during firing may cause damage and possible instrument malfunction.

CAUTION: Failure to advance the handle completely may result in incomplete staple formation, which may compromise the integrity of the staple line.

6. After firing, return firing knob all the way back to its pre-fire position.

7. Push in the cartridge release button to open the instrument (refer to step 1). Remove the instrument from the surgical site.

CAUTION: After firing, always inspect the staple line for hemostasis. Minor bleeding can be controlled by electrocautery or manual sutures.

NOTE: When a fired instrument is opened, the lockout device will deploy. The safety lockout prevents clamping of the instrument with a cartridge that has been fired.

RELOADING INSTRUCTIONS

1. To remove the fired cartridge, separate the instrument halves. Holding the cartridge half of the instrument, grasp the proximal end of the cartridge tabs and pull up and out to remove the cartridge from instrument.

2. To place a new reloading unit into the instrument, hold the reloading unit by the proximal end tabs and insert into the cartridge fork at a 30 to 45 degree angle from the distal end down until unit snaps into place. Remove shipping wedge after the reloading unit is fully loaded. This reloading unit will "float" up and down in final loaded position.

CAUTION: Ensure that a reloading unit is selected with the appropriate staple size for the tissue thickness. Overly thick or thin tissue may result in unacceptable staple formation.

STORE AT ROOM TEMPERATURE.

AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

DO NOT EXPOSE TO TEMPERATURES ABOVE 130° F (54° C).

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GOLD BRIDGE MEDICAL TECHNOLOGIES LTD
Shanghai International Holding Corporation GmbH